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UTILISATION CONFORME

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TB100 Control Solution Package Insert

INTENDED USE

The control solution is a desired aqueous glucose solution with small amounts of non-reactive ingredients and food coloring for verifying the accuracy of the entire blood glucose monitoring system. The control solution test examines the performance of the glucose meter, test strip and your operating techniques.

When to perform a control solution test

- When you are using your blood glucose monitoring system for the first time.
- When you are using a new batch of test strips.
- Anytime you question the performance of the system or on a regular basis to ensure accuracy, e.g. once a week.
- When you adjust your diabetic medication plans.
- When your blood glucose test result is lower or higher than your normal level.

Control solution contents

Each control solution is provided as a buffered, stabilized, aqueous solution of desired D(+) -glucose containing a preservative and other non-reactive ingredients. The solutions are formulated to react equivalently to a whole blood sample containing approximate 50, 110 or 300 mg/dL glucose concentrations. Actual lot-specific ranges are printed on the test strip bottle label.

Control Solution	mg/dL	mmol/L
Low	50	2,8
Normal	110	6,1
High	300	16,7

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
Tyson Bioresearch, Inc.
5F, No. 16,18,20,22, Kedong 3rd Rd., ZhuNan Township, Miaoli County 35053, Taiwan

Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41, D-30175 Hannover, Germany

Results falling outside the specified range may be caused by:

1. Error in performing the test
2. Control solution temperature is lower than 15°C (59°F) or higher than 35°C (95°F)
3. Expired or contaminated control solution
4. Expired or contaminated test strips
5. Meter malfunction

If the problem persists, the system may not be working properly. DO NOT continue to use the system to test your blood glucose level. If you continue to get results that fall outside the specified range, please contact Spengler Customer Service or e-mail: contact@spengler.fr

Storage and stability

For best results, store the control solution at room temperature; do not refrigerate. Refer to the expiry date for stability. The solution should NOT be used after the expiry date or 3 months after opening. For in vitro diagnostic use only.

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Testing

1. Insert a test strip, contact bar end first, into the test port. The meter will automatically turn on. A beep will sound and dotted lines moving clockwise will appear followed by a flashing .

2. After the  symbol appears, press the S button for 3 sec to switch to control solution mode. The  symbol indicates that the meter

Preparation

Make sure that you have set the glucometer in the control solution test "ctl" mode, so the result will not be stored in the memory and not be incorporated into the average. Please read your blood glucose monitoring system user manual.

Test procedure

Note: Make sure the solution is at room temperature (15°C to 35°C or 59°F to 95°F) prior to testing.

Comparison of control solution results

If the test result is outside the specified range printed on the test strip vial, repeat the test.

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SPENGLER SAS - 30 rue Jean de Guiramand