Holtex

# EASY 2

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Thank you for purchasing this Holtex sphygmomanometer. It's a practical tool that will let you check arterial blood pressure accurately.

### 1. What is blood pressure?

Arterial blood pressure is the pressure exerted by the blood against artery walls. Blood pressure changes continuously during the cardiac cycle. The highest level during the cycle is called systole. The lowest pressure is called diastole. Both of them have to be read to evaluate the condition of the arterial pressure. Many factors can influence blood pressure such as physical activity, anxiety or day time. Pressure is usually lower in the morning and becomes more elevated from afternoon until the evening. It's lower in summer than winter.

### 2. What is a normal arterial pressure?

Blood pressure varies from one person to another and depends on a number of factors such as age, weight and physical condition. The standard result for an adult between 18 and 45 years old is 120/180.

### 3. How is blood pressure measured?

#### - Installation in a suitable position

It's very important that the person is relaxed and comfortably seated. The arm must stay at the same height as the heart. The arm has to be extended without being tense. It should rest on a flat surface with the palm turned upwards. The arm shouldn't move during the measurement.

#### - How to fix the cuff?

Place the sphygmomanometer cuff two or three centimeters above the elbow on the bare arm. Fix it in position with Velcro parts. The cuff should be firmly held without being tight. One or two fingers can be introduced beneath the cuff. Never put it over clothes.

#### - Inflating the cuff

This product is provided with a push button. Press the button so as to remove all of the air from the cuff. Pump up the cuff using a regular rhythm until the manometer needle reaches 30 mmHG above the usual systolic pressure. If you don't know it, inflate to 200 mmHG. Introduce the head of a stethoscope beneath the cuff.

#### 4. Reading the systolic blood pressure

Press the valve slowly and keep the stethoscope head over the humeral artery. It is essential that an appropriate deflation speed is used to get precise reading. Do not keep the cuff inflated any longer than needed. When pumping starts, you should listen carefully with the stethoscope. Note the pressure shown on the manometer as soon as you start to hear weak rhythmic beats or muted beats. This is the systolic pressure reading.

#### 5. Reading the diastolic blood pressure

Let the pressure fall at the same speed as previously indicated. The diastolic pressure is reached when the muted beats stop. Deflate the cuff totally and remove stethoscope and it.

#### 6. Maintenance of the sphygmomanometer

If you clean it carefully, this product will give you satisfaction for many years. The cuff should be disinfected by being wiped with a suitable surface sanitizing solution.

Clean the gauge and bulb with a damp cloth.

Never inflate beyond 300 mm HG.

Never expose the cuff to strong sunlight.

Never touch the cuff with a sharp instrument.

Deflate the cuff before storing it.

Never disassemble the gauge.

Place the monitor in its storage pouch.

#### 7. Accuracy Control

Remove the tube from the manometer and hold the gauge upright. Your device is set correctly when the needle stops on the zero of the scale. If this is not the case, the instrument must be recalibrated. Contact your medical equipment reseller.

We recommend to check every 2 years the accuracy of your device.

## 8. Specifications

- Environmental conditions for which the minimum maximum tolerance of  $\pm 3$  mmHg is met:
  - Measurement conditions  $10^{\circ}$  to  $40^{\circ}$  C, relative humidity of 85%, without condensation
  - Storage Conditions  $-20^{\circ}$  to  $70^{\circ}$  C, relative humidity of 85%, without condensation
- Display Type: Round scale
- Scale graduation: increment of 2 mmHg
- Display range 0 to 300 mmHg
- Measuring range 0 to 280 mmHg
- Mobility of the needle: not reset
- Pressure generation: bulb
- Pressure reduction: air-release valve.

# EASY 2

30 rue Jean de Guiramand 13290 Aix en Provence  
Tél : +33(0)4 42 90 31 31 - Fax : +33(0)4 42 90 31 41  
[www.holtex.fr](http://www.holtex.fr) - [info@holtex.fr](mailto:info@holtex.fr)