# Spengler

## SPOT-CHECK 2 VS

**Station de diagnostic**

Diagnostic station

Estacion de diagnostico

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This manual is written for the current SPENGLER Spot-Check Monitor. The manual describes, in accordance with the Spot-Check Monitor's features and requirements, the main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc. as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for further details.

Version of This User Manual: Ver1

Issue Date: 6th November, 2018

Manufacturing date: See Label

All rights reserved.

Marks in the manual:

 **Caution: Instruction must be followed to avoid causing harm to the user or patient.**

 **Attention: must be followed to avoid causing damage to the Spot-Check monitor.**

 **Note: contains important information and advice about operations and application.**

## INSTRUCTIONS FOR USER

Dear Customers,

Thank you for purchasing the Spengler Spot-Check Monitor. Please read the following information before using the device. These instructions describe the operating procedures which are to be strictly followed, read these instructions carefully before using the Spengler Spot-Check Monitor. Failure to follow these instructions can cause monitoring abnormalities, damage to the monitor and personal injury. The manufacturer is NOT responsible for the safety, reliability and performance issues or any monitoring abnormalities, personal injury and equipment damage due to user's negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

### Warnings:

- **Do NOT use the device under flammable gas condition or in any environment that may lead to explosion**
- **The device and accessories that should not be serviced or maintained while the device is in use**
- **Do not modify this equipment without authorization from the manufacturer.**
- **The doctor or patient is the intended operator.**
- **The SpO<sub>2</sub> and NIBP measurements are frequently used functions.**
- **The device is IP22 and is protected against solid foreign objects of 12.5mm or greater, and protected against vertically falling water drops when the enclosure is tilted up to 15.**
- **Please check the monitor before use to verify that the accessories can function safely and correctly.**

- If the monitor is connected with other devices, the total leakage current may exceed the limitation and as a result this can cause potential danger
- Although biocompatibility tests have been performed on all the applied parts, under exceptional circumstances, allergic patients may have anaphylaxis. Do NOT use the monitor on patient with anaphylaxis
- All connecting cables and rubber tubes of the applied parts should be kept away from the patient's neck to prevent
- As a standard, please only use the components provided by the manufacturer or those that are of the same model and specifications as the accessories.
- If the monitor falls off a surface accidentally, please do NOT operate it before its safety and technical performance have been tested, and positive results obtained.
- Do NOT open the device cover without authorization. The cover should only be opened by a qualified service personnel.
- When disposing of the monitor and its accessories, the national regulation should be followed.
- There are some electromagnetic or inductance circuit designed in the device, use during MRI environment could burns or adversely affect the MRI image or the device's accuracy. So the device is MR unsafe.
- The device has no alarm and is intended only for spot-checking
- The device and accessories are provided non-sterile.

## TABLE OF CONTENTS

CHAPTER 1 OVERVIEW.....	43
1.1 Features.....	43
1.2 Product Name and Model.....	43
1.3 Intended Use.....	44
1.4 Impact on the Environment and Resources.....	44
CHAPTER 2 OPERATION INSTRUCTIONS.....	44
2.1 Appearance.....	44
2.1.1 The Front Panel.....	44
2.1.2 The Right and Upper Sides of the Device.....	46
2.2 Installation.....	47
2.2.1 Power Supply.....	47
2.2.2 Starting the Monitor.....	48
2.3 Taking Measurement.....	48
2.3.1 Blood Pressure Measurement.....	48
2.3.2 SpO <sub>2</sub> Measurement.....	52
2.4 Symbols.....	52
CHAPTER 3 MONITORING SCREEN DISPLAY.....	56
3.1 Measuring Screen.....	56
3.2 System Setting Screen.....	58
3.3 History Data Review Screen.....	60

CHAPTER 4 TECHNICAL SPECIFICATIONS.....	61
4.1 Blood Pressure Measurement.....	61
4.2 SpO <sub>2</sub> Measurement.....	62
4.3 Pulse Rate Measurement.....	63
4.4 Others.....	63
4.4.1 Operating Environment.....	63
4.4.2 Classification.....	64
CHAPTER 5 TROUBLESHOOTING.....	65
CHAPTER 6 ERROR MESSAGE INTERPRETING.....	66
CHAPTER 7 PACKING LIST.....	67
CHAPTER 8 MAINTENANCE AND SERVICE.....	68
8.1 Technical Maintenances.....	68
8.1.1 Daily Examination.....	68
8.1.2 Routine Maintenance.....	68
8.1.3 Battery Maintenance.....	68
8.2 Cleaning and Disinfection of the Main Unit.....	69
8.3 Cleaning and Disinfection of Accessories .....	70
8.4 Storage.....	71
8.5 Transportation.....	75

## CHAPTER 1 OVERVIEW

### 1.1 Features

- Small in size, light in weight, easy to carry and operate;
- Clear and large numeric display, segmented LCD panel, real-time clock display;
- Accurate blood pressure measurements can be activated or canceled by one shortcut button;
- Unique oximetry technique ensures quick and accurate SpO<sub>2</sub> & pulse rate measurements by smart sensors;
- Blood pressure, oxygen saturation, pulse rate and temperature can be measured simultaneously;
- Up to 100 user ID can be marked;
- Power management with power saving mode, auto power off and low battery indicator;
- Data storage with recall, up to 999 groups of records can be stored and recognised by patient ID.

### 1.2 Product Name and Model

Name: Spot-Check 2 VS

Model and Configuration:

Model	NIBP	SpO <sub>2</sub>	Pulse Rate	LCD display
PC-301	✓			✓
PC-302	✓	✓	✓	✓

### 1.3 Intended Use

The Spot-Check Monitor is a device designed for spot-checking the user's physiological parameters, such as non-invasive blood pressure (NIBP), functional oxygen saturation (SpO<sub>2</sub>) and pulse rate (PR). This device is applicable for use in clinical institutions and has no conditions or factors of contraindication.

### 1.4 Impact on the Environment and Resources

Low

## CHAPTER 2 OPERATION INSTRUCTIONS

### 2.1 Appearance

#### 2.1.1 The Front Panel


Figure 2.1

**Description :**

1/2.  Up/down key: on the setup display screen, a short press will change the parameter value step by step, press and hold to change the parameter values quickly; on the review display screen, short press to review the history data records, press and hold to recall the history data records quickly.

3.  Memory key: on the measurement display screen, press and hold the key (for 3 seconds) to enter into the review display screen; once the review display screen, a short press will recall the history data records.

On the setup display screen, all parameters can be set in anticlockwise order by pressing and holding  key, similarly, a short press of  key will set the parameters in clockwise order.

4.  Menu key: on the measurement display screen, press and hold the menu key to enter the setup screen; on the setup or review display screen, press and hold the  key to go back to the measurement display screen.

5.  Start/cancel button: on the measurement display screen, a short press of this button will activate or cancel the blood pressure measurement.

### 2.1.2 The Right and Upper Sides of the Device


Figure 2.2


Figure 2.3

The power switch and external DC power input socket are on the right side of the monitor as shown in figure 2.2

The signal input/output ports are on the upper side of the monitor as shown in figure 2.3.

**Description :**

DC 5.0V 1.2A

1.  : External DC power input socket.

2.  : Power switch: = press and hold
3. Battery cover.
4.  : Charge / USB data interface.
5. Unused
6. NIBP: Cuff connector.
7. SpO2: Connector to link with the smart SpO2 probe.
8. Unused

## 2.2 Installation

### 2.2.1 Power Supply

1. Internal power supply to the built-in battery

When the battery indicator  displays full grids, the built-in battery is fully charged. When it blinks, the battery voltage is low, and the user should charge the battery by connecting the device to the AC power adapter. When the grids of the battery indicator are rolling circularly, the battery is being charged.

2. External power supply from the AC power adapter Use the AC power adapter provided by the manufacturer. Make sure that the mains power supply is 110~240VAC with 50/60Hz.

## 2.2.2 Starting the Monitor

By pressing and holding down the switch, the software version will be displayed after releasing the switch, the device will enter the measurement display screen automatically. The user can then begin to operate the monitor.

Note: Do not touch the 4 touch-keys during switching on the monitor, or the touch-keys may be out of work temporarily. If the touch-keys are not sensitive enough, do not operate them for over 8 seconds, then the touch-keys will resume to its normal sensitivity automatically.

- The monitor is powered by built-in Lithium battery, when the battery voltage is low, the measurement and wireless connection may be unstable.
- For electric safety, do not use the monitor during its battery is in charging.
- Please make measurement for a single person at a time.

> If monitor fails to start by pressing the switch, please use the external power supply.

## 2.3 Taking Measurement

### 2.3.1 Blood Pressure Measurement

1. An appropriate cuff should be selected according to the age and arm circumference of the patient. The cuff width should be 2/3 of the length of the upper arm. The inflatable part should be long enough to permit wrapping approximately 80% of the limb.
2. Applying the cuff: unfold the cuff and wrap it around the upper arm evenly to the appropriate tightness. The correct cuff position is shown in figure 2.5



Figure 2.5 Cuff position

3.3. Connect the hose from the cuff to the connector on the upper-side of the device where marked "NIBP".

4. Long press menu key «» to enter the setup screen, select the proper patient category, for example, if you select «Child», then the indicator on the screen will focus towards «Child»,

shown as «», which means that the current patient type is set to Child. Note: refer to Section 3.2 for detailed setting operation.

5. Press the start/cancel button «» to begin the blood pressure measurement. The result will be displayed on the screen, and the corresponding blood pressure level will be indicated,

shown as «», which means that the adult's blood pressure result is normal. Please refer to the form below for detailed blood pressure levels. Note: blood pressure level indication is only for «Adult» patient type.

<b>Level</b>	<b>Blood pressure range (Unit: mmHg)</b>
N (Normal)	SYS<130mmHg, and DIA<85mmHg ;
H-N (High normal)	130mmHg≤SYS<140mmHg, and 85mmHg ≤ DIA<90mmHg ;
HT (Hypertension)	SYS≥140mmHg, or DIA≥90mmHg ;

### **NIBP Measurement Principle**

The NIBP measurement is based on oscillation technology. The measurement is started by User Manual for Spot-Check Monitor 14 inflating the cuff by a pump automatically after the cuff pressure is high enough to block the arterial blood flow within the upper arm, then the cuff pressure is deflated slowly, and all the change of cuff pressure in the deflation process is recorded to calculate blood pressure based on certain algorithm. The device will judge whether the quality of signal is good enough. If the signal is not good enough (such as sudden movement or touch of cuff while measurement), the device will stop deflating or re-inflating, or aborting this measurement and calculation.

### **Safety Instructions for blood pressure measurement**

- Blood pressure measurement is prohibited to those who have severe hemorrhagic tendencies or with sickle cell disease, as partial bleeding maybe caused.
- Too frequent measurements or connection tube kinking may result in purpura, neuralgia and lack of blood.
- Wrap the cuff and operation of the start/cancel button  are the frequently used functions.
- Do NOT apply the CUFF over a wound, as this can cause further injury.
- Operation of the device does not result in prolonged impairment of PATIENT blood

circulation.


- Do NOT wrap the cuff on limbs with transfusion tubes, intubation or skin lesions on the area, as damage may be.
- The equipment can be used on pregnant or pre-eclamptic patients.
- The proper operating steps are needed to obtain accurate resting Blood Pressure reading routinely:
  - Patient position in normal state, including comfortably seated, legs uncrossed, feet flat on the floor, back and arm supported, middle of the cuff at the level of the right atrium of the heart.
  - The patient should be as relaxed as possible and should not talk during the measurement procedure.
  - 5 minutes should elapse before the first reading is taken.
- The user needs to check the operation of this equipment does not result in prolonged impairment of patient blood circulation.
- Readings can be affected by the measurement site, the position of the patient(standing, sitting, lying down), exercise, or the patient's physiological.
- The performance of the equipment can be affected by extreme temperature, humidity and altitude.
- Avoid compressing or restricting the connection tubing.
- The patient should be comfortably seated with their legs uncrossed and feet flat on the floor. The patient's back and arm should be supported, and the middle of the cuff should be level with the right atrium of the heart.
- The environment or operational factors which can affect the performance of the device and/or its blood pressure reading (e.g. common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, arterial sclerosis, poor perfusion diabetes, age, pregnancy, pre-eclampsia, renal diseases, patient motion, trembling, shivering).

- ⚠ Measurements should be taken at appropriate intervals. Frequent measurements with short intervals may lead to pressed arm, reduced blood flow, low blood pressure, and result in an inaccurate reading. It is recommended that the measurements are taken in intervals of more than two minutes
- ⚠ Before use, empty the cuff until there is no residual air inside.
- ⚠ Do NOT allow the cuff to twist or bend.
- ⚠ Do NOT twist the cuff hose or put heavy things on it.
- ⚠ Please hold the connector of the hose while connecting and disconnecting it to the device.
- ⚠ If arrhythmia or auricular fibrillation occurs, take measurement again.
- ⚠ The patient should sit or lay down with calm condition and make the cuff and the patient's heart on the same level to get accurate measurement. Other positions may lead to inaccurate measurement.

### 2.3.2 SpO<sub>2</sub> Measurement

#### Operation procedures:

1. Connect the smart SpO<sub>2</sub> probe to the connector on the upper-side of the device marked «SpO<sub>2</sub>». When disconnecting the connector, be sure to hold the head of the connector firmly and pull
2. The red blinking light inside the clip of the SpO<sub>2</sub> probe indicates a successful connection.
3. Insert one finger (index finger is preferred, the nail should be not too long) into the clip of the probe according to the finger mark shown as below.
4. The device will begin to take the measurement.

### **SpO<sub>2</sub> Measurement Principle**


SpO<sub>2</sub> measurement is based on dual wavelength opto-plethysmometry technology, a unique hardware and software design as its proprietary technology from Shenzhen Creative Industry Co., Ltd.. By use of red and infra-red light emitting through the patient's finger, the photo-detector at the other side senses the transmitted light and converts to current for later amplification and filtering. The acquired light intensity signals (plethysmogram) are digitalized and further processed with proprietary algorithm to determine the SpO<sub>2</sub> and pulse rate value.

### **Safety instructions for SpO<sub>2</sub> measurements**

- ⚠ Continuous use of the SpO<sub>2</sub> probe may result in discomfort or pain, especially for those with microcirculatory problems. It is recommended that the probe should NOT be applied to the same place for over two hours, change the measurement site periodically and when necessary.
- ⚠ Do NOT place the SpO<sub>2</sub> probe on a finger with edema or fragile tissue.

- ⚠ When the ambient temperature is over 35°C, please change the measuring site every two hours; when the ambient temperature is over 37°C, please do NOT use the SpO2 sensor, as using in high temperatures can cause burns
- ⚠ Do NOT put the SpO2 probe and pressure cuff on the same limb, otherwise the blood pressure measurement may affect the SpO2 measurement.
- ⚠ The device is calibrated to display functional oxygen saturation.
- ⚠ Do NOT allow the sensor cable to twist or bend.
- ⚠ Check the SpO2 sensor and cable before use. Do NOT use a damaged SpO2 sensor
- ⚠ When the temperature of the SpO2 sensor is abnormal, do not use it further
- ⚠ Remove nail polisher or other cosmetic products from the fingernail.
- ⚠ The fingernail should be of normal length.
- ⚠ The SpO2 sensor cannot be immersed into water, liquid or cleanser.
- ⚠ The SpO2 sensor can be repeatedly used. Please clean and disinfect before reuse.
- ⚠ The SpO2 sensor can be repeatedly used. Please clean and disinfect before reuse
- ⚠ Anemia or low hemoglobin concentrations, intravascular dyes, carboxyhemoglobin, methemoglobin, and dysfunctional hemoglobin may effect the SpO2 accuracy. If the patient has such situation, do not rely on the measured result for diagnostic decision, and it's recommended for the patient to consult with the doctor.

## 2.4 Symbols

Symbol	Description	Symbol	Description
	Wireless		Pulse rate (unit: bpm, beat per min)
	Alarm		Battery voltage indicator
	Memory icon		USB icon
	Power on/off switch		External DC power input
	Type BF applied parts		Charger or USB Data interface / Connector to link with blood glucose meter
	Serial Number		Refer to manual
	Manufacturer information		No SpO2 Alarms
	Class II device		Follow WEEE regulations for disposal
	Batch Number		

## CHAPTER 3 MONITORING SCREEN DISPLAY

### 3.1 Measuring Screen

#### Screen Description:

1. : USB connection icon.

2. : Non available in PC 301 et PC 302 versions.

Wireless transmission icon means that the wireless transmission function is on; when the icon is blinking, the wireless connection set up is unsuccessful; when this icon is steady, the wireless connection set up is successful ; : the wireless transmission function is off.

3. : Beep sound indicator; : pulse beep is on; : pulse beep is off.
4. : Battery voltage indicator. When the battery is full, the battery voltage indicator displays a full grid. When the indicator is blinking, it means the battery voltage is low and the user should charge the battery. Please connect the device to the external power supply to ensure the correct use of the monitor, and the battery will be charged. During charging, the grids in the battery indicator will roll circularly.
- 5~7. : Means the inflation pressure during cuff inflation. When displaying the measurement result, the indicator appears on the corresponding blood pressure level, such as N (Normal), H-N (High normal), and HT (hypertension).
- 8~10 : indicates the NIBP measurement mode, the indicator appears on the corresponding patient type, for Adult, for Child, and for Infant.
11. M: Memory
12. ID: the patient ID, which can be set from 0 to 99.
13. NO.: the number of stored data, up to 999 records can be stored for each ID.
14. H:M: the time stamp (hour: minute). The time can be set in the system setup screen
15. M-D: the time stamp (month-day). The date can be set in the system setup screen.
16. SYST: Systolic pressure
17. DIAS: Diastolic pressure
18. kPa/mmHg: unit of blood pressure, 1kPa=7.5mmHg.
19. SpO2: the value of SpO2 with unit of %.
20. : Pulse bar-graph.
21. PR: Pulse rate with unit of bpm.
22. : the heart beat symbol which flashes with heart beat.

### 3.2 System Setting Screen

On the measurement display screen, press and hold the menu key to setup the display screen, as shown in figure 3.2. The user can choose the settings for the wireless function, pulse beep, blood pressure unit, , date and time.

**Operation Instruction:**

1. Press and hold the «» key, and release after hearing one beep, to enter into the setup screen. When the patient ID blinks, the setup function is available.
2. A short /long press of the «» key to cycle select the setting item, the item will be blinking if selected. For example, if NIBP patient type is selected, then NIBP indicator will be blinking.  
(All parameters can be set in the order of anti-clockwise by pressing and holding the «» key, the order will be: ID -> Minute -> Hour -> Day -> Month -> Year -> TEMP -> NIBP unit -> Alarm sound -> Wireless transmission function -> NIBP mode).
3. A short press of the «/» key sets the detailed item. For example, if you want to choose Child type, then press up/down key to move the indicator on the label of «Child».
4. A short press of the «» key confirms the setting.
5. Press and hold the «» key brings the screen display back to the measurement display screen. The monitor will also switch back to the measurement display screen as well if there has been no operation for 30 seconds.

Note:

1. On the setup display screen, all parameters can be set in anticlockwise order by pressing and holding the «» key.
2. For setting the date, the century is fixed to be 20, i.e. "13y" indicates the year 2013. Please see the following example for the date and time: 11:14», March 23, 2013.

13  
11:14  
03-23

### 3.3 History Data Review Screen

On the measurement display screen, press and hold «» key to recall the stored data records, as shown in figure 3.

Figure 3

**Operation instructions:**

1. Press and hold the «» key, release the key after hearing one beep. The memory mark «» will appear (i.e. entering to review display screen). The patient's ID number will blink at the same time.
2. A short press of the «/» key will browse the patient's ID number.
3. A short press of the «» key will confirm the setting, and the recorded number (No.). will blink.
4. A short press of the «/» key will set the recorded number to be recalled. The data displayed on the screen is for the specific record of the patient that you selected.

Note: when selecting the patient ID, the screen only displays patients with history data records.

## **CHAPTER 4 TECHNICAL SPECIFICATIONS**

### **4.1 Blood Pressure Measurement**

1. Technique: Oscillometric
2. Pressure measuring range: 0mmHg~300mmHg
3. Overpressure protection limit  
Adult: 300mmHg(39.9kPa); Child: 240mmHg; Infant: 150mmHg
4. Blood pressure measurement range:  
SYST: 30mmHg~270mmHg DIAS: 20mmHg~235mmHg MAP: 10mmHg~220mmHg
5. Blood pressure measuring accuracy:  
Maximal mean difference:  $\leq \pm 5\text{mmHg}$  (0.67kPa)  
Maximal standard deviation:  $\leq 8\text{mmHg}$  (1.067kPa)

## 4.2 SpO<sub>2</sub> Measurement

1. Technique: optical with dual-wavelength LED wavelength: Red light: 663nm, Infrared light: 890nm  
Maximal optical output power: less than 2mW maximum average
2. SpO<sub>2</sub> display range: 0%~100%
3. SpO<sub>2</sub> measuring accuracy: Arms is not greater than 3% for SpO<sub>2</sub> range from 70% to 100%, undefined during range 0%~70%  
Note: Arms is defined as root-mean-square value of deviation according to ISO 80601-2-61 / ISO 9919
4. Measuring mode: spot-checking
5. SpO<sub>2</sub> display update: every second
6. SpO<sub>2</sub> averaging: Averages the recent eight seconds value falling within the acceptable limits
7. The table with measured SpO<sub>2</sub> accuracy specification in the discrete SpO<sub>2</sub> ranges:

SpO <sub>2</sub> range	Arms
70%~80%	1.62
80%~90%	1.09
90%~100%	1.58
70%~100%	1.42

### 4.3 Pulse Rate Measurement

1. PR measuring range: 30bpm~240bpm
2. Pulse rate measuring accuracy:  $\pm 2\text{bpm}$  or  $\pm 2\%$ , which is greater

### 4.4 Others

#### 4.4.1 Operating Environment

1. Operating temperature: 5°C~40°C; Relative humidity:30%~80%; Atmospheric pressure: 70.0kPa~106.0kPa; Power supply: a.c. 110V-240V AC, 50/60Hz;  
Internal power supply: d.c.3.7V(rechargeable Lithium battery); Input: 15VA
2. The device should be situated in a place protected against direct sunlight, to prevent overheating inside of the equipment.
3. Do not use this equipment in combination with any equipment other than those approved in the user manual.
4. The device should be stored and used in a specified temperature, humidity and atmospheric pressure range, or damage maybe caused to the device
5. If the device gets wet by accident, the operator should NOT turn on the power until it has been thoroughly air dried.
6. Do not use this equipment in an environment with toxic or inflammable gas.
7. Only monitor a single person at a time.
8. Do not expose the device to a magnetic resonance (MR) environment
  - The device may present a risk of projectile injury due to the presence of ferromagnetic materials that can be attracted by the MR magnet core.
  - Thermal injury and burns may occur due to the metal components of the device that can heat during MR scanning

- The device may generate artifacts in the MR image.
- The device may not function properly due to the strong magnetic and radiofrequency fields generated by the MR scanner.

Warning: Do not use any other adapters than those provided by SPENGLER.

#### **4.4.2 Classification**

1. Protection against electric shock: Class II equipment and internally powered equipment
2. The degree of protection against electric shock: Type BF applied part
3. Define apply part: cuff, SpO<sub>2</sub> probe.
4. The degree of protection against harmful ingress of liquid: The equipment is IP22 with protection against ingress of liquid 5. Electro-magnetic Compatibility: Group I, Class A

## CHAPTER 5 TROUBLESHOOTING

Trouble	Possible reason	Solution
Cannot turn on the device	The built-in battery is drained	Recharge by connecting the power supply adapter
	Battery is not installed	Install the lithium battery
	Some parts provided by other manufacturers are inserted to the connector	Remove the related parts and try again
No blood pressure result	The cuff is wrapped around the arm incorrectly	Wrap the cuff around the arm correctly
	The windpipe is not correctly inserted to the NIBP jack	Insert the windpipe to the NIBP jack
No SpO2 result	The SpO2 probe is not plugged into the «SpO2» connector (or Port1 or Port2)	Plug the SpO2 probe into the «SpO2» connector (or Port1 or Port2)

## CHAPTER 6 ERROR MESSAGE INTERPRETING

Error code	Description
ERR 01	Fails to inflate pressure to 30 mmHg within 7 seconds (The cuff is not well-wrapped)
ERR 02	Cuff pressure is over 295 mmHg (Overpressure protection)
ERR 03	No valid pulse is detected
ERR 04	Excessive motion artifact
ERR 05	Invalid measured result
ERR 06	Air leakage is detected
ERR 07	Self-test fail
ERR 08	Pneumatic system error
ERR 09	Signal saturation
ERR 10	Air leakage check fail
ERR 11	Hardware fault
ERR 12	Measurement timeout (the max. measurement time for adult is 120 seconds when the pressure is over 200 mmHg, 90 seconds when the pressure is less than 200 mmHg, and 90 seconds for neonate.)

## CHAPTER 7 PACKING LIST

Item	Description	Quantity	Check
1	Spot check monitor	1	OK
2	Handbag	1	OK
3	User manual	1	OK
4	Cuff	1	OK
5	Charger	1	OK
6	SPO2 probe	1	OK

## CHAPTER 8 MAINTENANCE AND SERVICE

The Spot-Check Monitor should be properly maintained to ensure its maximum performance and long service life. In addition to the warranty period, the company also offers long-term service for each customer. It is important that the user reads and follows the operation instructions, important information and maintenance measures.

### 8.1 Technical Maintenances

**8.1.1 Daily Examination :** Before using the monitor, the following checks should be carried out:

- Check the monitor for any mechanical damage;
- Inspect the exposed parts and the inserted parts of all the leads, and the accessories;

- Examine all the functions of the monitor that are likely to be used for patient monitoring, and ensure that it is in good working condition. If there is any indication of damage, or if damage is accurately proven, do not use the device. Contact your supplier for advice and to reach a satisfaction solution.

### **8.1.2 Routine Maintenance**

- ⚠ If the hospital fails to carry out a satisfactory maintenance program on the monitor, it may cause harm to the patient.
- ⚠ If there is any indication of cable and transducer deterioration or damage, please do not use.
- ⚠ The SpO<sub>2</sub> function has been adjusted before vending. If the user needs to adjust the SpO<sub>2</sub>, adjust by using the simulator mode FLUKE INDEX2.
- ⚠ The adjustable units in the monitor such as the potentiometer cannot be adjusted without permission so to avoid unnecessary failures that may affect the normal application.
- ⚠ It is recommended to use the battery once a month to ensure a strong power supply capacity and long service life, and recharge once the power has completely run out.

### **8.1.3 Battery Maintenance**

- ⚠ Please pay attention to the polarity of the battery, do NOT insert into the battery compartment with reversed polarities.
- ⚠ In order to avoid damaging the battery, do NOT use other power supply device to charge the battery
- ⚠ After use, dispose of the battery according to local regulations, do NOT throw into fire.
- ⚠ Do NOT hit or strike the battery with force

- ⚠ Do NOT use this battery in other devices
- ⚠ Do NOT use this battery below -20°C or above 60°C.
- ⚠ In order to maintain the battery supply and prolong the battery lifetime, please charge the battery routinely. Regularly, charge the battery every 3 months even if the device has not been used.
- ⚠ Only use a battery with the specification recommended by the manufacturer.
- ⚠ Whether the monitor is on or off, the built-in battery will charge as long as the monitor is connected to an AC adapter and the AC power is on. When the battery is full, it will stop charging to avoid causing any damage. If the monitor is connected to an AC adapter and the AC power is on, it will use the AC power, but when the AC power is off, the battery power will be used. Priority of using the AC power and the power switch between the AC and battery is automatic and seamless.
- ⚠ If the battery is damaged, please replace it with a battery with «CCC» or «CE» mark. The model and specifications of the battery should be the same as the original battery. The user must ensure that the battery meets all applicable safety codes. The user can also contact the distributor for service

## 8.2 Cleaning and Disinfection of the Main Unit

- ⚠ Switch off the monitor and disconnect the power cord before cleaning.
- ⚠ Keep the monitor free from dust
- ⚠ It is recommended to regularly clean the outer shell and screen of the monitor. Only use a non-corrosive cleanser such as clear water.

- ⌚ Wipe the surface of the monitor and transducers with an alcohol impregnated wipe, and dry with a clean cloth or just air-dry.
- ⌚ Dilute the cleaner.
- ⌚ Do NOT use scrubbing materials.
- ⌚ The monitor can be disinfected. Please clear the monitor first
- ⌚ To avoid damage do not let liquid cleaner flow into the connector jack of the monitor
- ⌚ Clean the exterior of the connector only.
- ⌚ Do NOT let any liquid flow into the shell or any other parts of the monitor
- ⌚ Do NOT leave any residue liquid or disinfectant on the surface of the monitor.
- ⌚ Do NOT perform high pressure sterilization on the monitor
- ⌚ Do NOT immerse any parts of the monitor or its accessories in the liquid.
- ⌚ If the monitor accidentally becomes wet, it should be thoroughly dried before use. The rear cover can be removed by a qualified service technician to verify the absence of water
- ⌚ Do NOT pour disinfectant on the monitor's surface while disinfecting

### **8.3 Cleaning and Disinfection of Accessories**

It is recommended to clean and disinfect the accessories (excluding the SpO<sub>2</sub> probe) with a piece of gauze soaked in 75% Alcohol or 70% Isopropanol.

- ⌚ Do not use damaged accessories.
- ⌚ Accessories cannot be entirely immersed into water, liquid or cleanser.
- ⌚ Do NOT use radiation, steam or epoxyethane to disinfect accessories.

- ⚠ Wipe off any remaining residue of alcohol or isopropanol after disinfection.
- ⚠ To prevent cross infection, the user wipes the temperature sensitive probe with 75% alcohol before and after use, then wipes the residue clean with clean dry cloth

#### **8.4 Storage**

If the equipment will not be used for a long time period of time, wipe it clean and return it to the packaging. Store in a dry well ventilated place free from dust and corrosive gases.

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Pression atmosphérique : 53,0 kPa~106,0 kPa

#### **8.5 Transportation**

The monitor should be transported by land (vehicle or railway) or air in accordance with the contractual terms. Do NOT hit or drop with force.

#### **Appendix I Classification of Blood Pressure Level**

**Classification of Blood Pressure Level**

Category	SYS (mmHg)	DIA (mmHg)
Optimal	<120	<80
Normal	<130	<85
High normal	130~139	85~89
Grade 1 hypertension (Mild hypertension)	140~159	90~99
Grade 2 hypertension (Moderate hypertension)	160~179	100~109
Grade 3 hypertension (Severe hypertension)	≥180	≥110
Isolated systolic hypertension	≥140	<90

Reference : The 1999 WHO-ISH Guidelines for the Management of Hypertension

## **Warranty Clause**

1. We offer one year free repair for the main unit and 6-month for the accessories since purchasing. The invoice date is the starting date.
2. Make sure that the Warranty is sealed. Please present both the Warranty and the Invoice when you are Disassemble or refit the device arbitrarily; requesting free repair.
3. Please send the device to the specified place for repair.
4. We will not offer free repair if one of the following condition exists: Incorrect operation; Drop, immerse and collision of the device; Incorrect repair by unauthorized personnel.
5. The cost of repair beyond warranty will be taken by the charge rate, and we offer 3-month repair for such malfunction.

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Patent State Intellectual Property Office of the P.R.C. patented and certificated Creative Spot-Check Monitor on March 26 th , 2014.

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