# Spengler

## LIAN® DESSILLONS

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Sphygmomanometer/Tensiómetro aneroide

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|                                                        |    |
|--------------------------------------------------------|----|
| 1. PRODUCT IDENTIFICATION .....                        | 12 |
| 2. INTENDED USE .....                                  | 14 |
| 3. CONTRAINDICATIONS, COMPLICATIONS AND CAUTIONS ..... | 15 |
| 4. TECHNICAL SPECIFICATIONS .....                      | 15 |
| 5. CALIBRATION AND VERIFICATION .....                  | 16 |
| 6. USE AND PERFORMANCES OF THE DEVICE .....            | 17 |
| 7. CLEANING, STORAGE AND DISPOSAL .....                | 18 |
| 8. INFORMATION IN CASE OF SERIOUS INCIDENT .....       | 19 |
| 9. WARRANTY .....                                      | 19 |
| 10. TABLE OF USED SYMBOLS .....                        | 19 |

## 1. PRODUCT IDENTIFICATION

The sphygmomanometer is a medical device that measures the pressure exercised by blood on arteries by the auscultatory method requiring a cuff or lower limb cuff and a stethoscope.

The sphygmomanometer is a mano-bulb type with aneroid diaphragm, very light.


The precision pressure gauge graduated from 0 to 300mmHg is fixed on an anodized aluminium valve, equipped with a decompression knob.

The readings of measured valued is facilitated by the size of the dial. Guaranteed pressure accuracy of  $\pm 3$ mm over the whole measuring range.

- The medical class of the device is Im
- Lifetime is subjected to device use and cleaning and decontamination cycles. Usual lifetime is 3 years for the sphygmomanometer.
- This device is intended to be used with other devices such as: cuffs, lower limb cuffs (Soft Cuff® and CLINIC® brands) and stethoscopes to determine arterial pressure.
- The size of the cuff has to be adapted to patient's morphology.


### LIAN® DESSILLONS references

| DESIGNATION                                                         | REFERENCE | WITH Soft Cuff® | WITH CLINIC® CUFF |
|---------------------------------------------------------------------|-----------|-----------------|-------------------|
| Sphygmomanometer Lian® Dessillons Soft Cuff® small adult SPENGLER   |           | TS50003S        |                   |
| Sphygmomanometer Lian® Dessillons Soft Cuff® adult SPENGLER         |           | TS50004S        |                   |
| Sphygmomanometer Lian® Dessillons Soft Cuff® large adult SPENGLER   |           | TS50005S        |                   |
| Sphygmomanometer Lian® Dessillons Soft Cuff® multi cuffs SPENGLER   |           | TS53000S        |                   |
| Sphygmomanometer Lian® Dessillons CLINIC® small adult cuff SPENGLER |           |                 | TI50003S          |
| Sphygmomanometer Lian® Dessillons CLINIC® adult cuff SPENGLER       |           |                 | TI50004S          |
| Sphygmomanometer Lian® Dessillons CLINIC® large adult cuff SPENGLER |           |                 | TI50005S          |
| Sphygmomanometer Lian® Dessillons CLINIC® multi cuffs SPENGLER      |           |                 | TI53000S          |
| Sphygmomanometer Lian® Dessillons single manometer SPENGLER         | A10370S   |                 |                   |


## References and sizes of the Soft Cuff® cuffs and lower limb cuff

| REFERENCE WITH CONNECTOR | DESIGNATION                 | OUTPUT | CIRCUMFERENCE (CM) |      | POCKET SIZE (CM) |
|--------------------------|-----------------------------|--------|--------------------|------|------------------|
|                          |                             |        | MIN                | MAX. |                  |
| BT101S                   | NIBP Soft Cuff® newborn     | 1      | 6                  | 12   | 4                |
| BT102S                   | NIBP Soft Cuff® child       | 1      | 11                 | 22   | 8                |
| BT103S                   | NIBP Soft Cuff® small adult | 1      | 16                 | 28   | 11               |
| BT104S                   | NIBP Soft Cuff® adult       | 1      | 26                 | 33   | 13               |
| BT105S                   | NIBP Soft Cuff® large adult | 1      | 33                 | 41   | 15               |
| CT104S                   | NIBP thigh Soft Cuff® adult | 1      | 39                 | 55   | 19               |

## References and sizes of the CLINIC® cuffs and lower limb cuff

| REFERENCE WITH CONNECTOR AND CAP | DESIGNATION                              | OUTPUT | CIRCUMFERENCE (CM) |      | POCKET SIZE (CM) |
|----------------------------------|------------------------------------------|--------|--------------------|------|------------------|
|                                  |                                          |        | MIN                | MAX. |                  |
| BM101S                           | NIBP cuff immersible CLINIC® newborn     | 1      | 6                  | 12   | 4                |
| BM102S                           | NIBP cuff immersible CLINIC® child       | 1      | 11                 | 22   | 8                |
| BM103S                           | NIBP cuff immersible CLINIC® small adult | 1      | 16                 | 28   | 11               |
| BM104S                           | NIBP cuff immersible CLINIC® adult       | 1      | 26                 | 33   | 13               |
| BM105S                           | NIBP cuff immersible CLINIC® large adult | 1      | 33                 | 41   | 15               |
| CM104S                           | NIBP thigh cuff immersible CLINIC®       | 1      | 39                 | 55   | 19               |

## 2. INTENDED USE

|                            |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>INTENDED USE</b>        | <p>Sphygmomanometers are designed for measuring non invasive blood pressure. The sphygmomanometer allows to measure the pressures (systolic and diastolic) of blood circulating in arteries.</p> <p>The readings are expressed in mmHg (millimeters of mercury) and include 2 values:</p> <ul style="list-style-type: none"><li>• Systolic pressure</li><li>• Diastolic pressure</li></ul> <p>The blood pressure measurement is done with an adapted cuff to patient's morphology and a stethoscope. It then allows a medical diagnosis if needed.</p> |
| <b>INDICATION</b>          | <p>Medical diagnosis, post-surgical surveillance, pregnant woman, patients with chronical disease (hypertension, diabetes , etc...).</p> <p>The non invasive measurement of blood pressure is done with a sphygmomanometer and a stethoscope. The sphygmomanometer allows to measure the back pressure on an artery coming from the inflated cuff while examining the artery.</p>                                                                                                                                                                      |
| <b>PATIENTS POPULATION</b> | <p>Newborns, pediatric patients, adults and elderly people.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>USERS</b>               | <p>Users are healthcare professionals trained to the non invasive measurement of blood pressure technique using the auscultatory method.</p> <p>The device is not intended fo use by profane people.</p>                                                                                                                                                                                                                                                                                                                                               |

### 3. CONTRAINDICATIONS, COMPLICATIONS AND CAUTIONS

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>CONTRAINDICATIONS</b> | No known contraindication                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>COMPLICATIONS</b>     | Side effect linked to the blood pressure measurement: none                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>CAUTIONS</b>          | <p><b><u>Cautions before use:</u></b></p> <ol style="list-style-type: none"> <li>1. Select the cuff adapted to patient's morphology</li> <li>2. It is an absolute must that the cuff features the manometer adapted connector.</li> <li>3. A check of the material is necessary before use: visually control the overall condition, and particularly the presence of the manometer, the decompression valve, the inflation bulb. If the device shows a dysfunction (cuff or manometer leakage), it has not to be used;</li> </ol> <p><b><u>Cautions while using and after use:</u></b></p> <ol style="list-style-type: none"> <li>4. In case of device leakage, inflate again the cuff or change the whole device.</li> <li>5. Cuff and lower limb cuff should not be manipulated with any sharp objects</li> <li>6. Comply with cleaning and disinfection instructions, check the device after disinfection.</li> </ol> |

### 4. TECHNICAL SPECIFICATIONS

|                        | <b>MATERIAL</b>                                                                                                                      | <b>PRESSURE READING</b>                                            | <b>FASTENING SYSTEM</b>                |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|----------------------------------------|
| <b>Manometer</b>       | Aluminium body connector<br>PVC bulb<br>Polycarbonate glass<br>Manometer protection: TPU<br>Decompression valve: nickel-plated brass | Graduated manometer<br>0 - 300 mmHg,<br><br>Precision $\pm 3$ mmHg | Female (CMP001137) on manometer        |
| <b>Soft Cuff® Cuff</b> | PU leather<br>Bladder in welded PU                                                                                                   | N/A                                                                | Tube with A11313 (cuff male connector) |
| <b>CLINIC® Cuff</b>    | PU coated fabric, high frequency welded                                                                                              | N/A                                                                | Tube with A11313 (cuff male connector) |

## 5. CALIBRATION AND VERIFICATION

- Refer to the local regulation for calibration frequency of the user. The calibration must be carried out by SPENGLER MEDICAL's customer service or an authorized person by SPENGLER MEDICAL. To request technical support for SPENGLER MEDICAL, call +33 (0)4 42 90 31 31 or send an email to SAV@spengler-med.fr.

INSTRUCTIONS  
FOR USE

1. It is important to check that the chosen Spengler cuff has the appropriate connector for the tensiometer.
2. Select the cuff adapted to patient's morphology Size and usage range are indicated on the cuff.
3. Leave the patient rest for 5 minutes prior to test.
4. Place the patient in a comfortable position, back and arm supported, feet uncrossed and flat on the floor without moving nor speaking.
5. Position the center of the bladder represented by a heart on the axis of the humeral artery.
6. Wrap the cuff properly around the arm and ensure its retention using the intended fastening system. The cuff should be positioned on the upper arm at the level of the right atrium.
7. Once the cuff is positioned and ready for use, the INDEX must be between the Min and the Max (heart axis). If the INDEX exceeds the Max index at the heart axis, replace the cuff with a larger size model. If the INDEX exceeds the Min index at the tip of the range, replace the cuff with a smaller size model.
8. When the cuff and INDEX are correctly positioned the blood pressure measurement process can begin.
9. Position the tips of the stethoscope lyre to your ears, the chestpiece in the listening position, just below the cuff, on the path of the artery.
10. Ensure that the pear pressure relief screw is in the closed position and inflate to 30 or 40 mmHg above the assumed pressure. The blood no longer circulates at the microphone, no sound is audible at the stethoscope.
11. Slightly unscrew the bulb decompression screw to slowly reduce the inflation pressure and allow blood to flow again, the first beats are noticeable from the stethoscope.
12. Read at this moment the pressure indicated on the manometer, this value corresponds to the maximum pressure or systolic pressure.
13. Continue decompression, the beats perceived with the stethoscope amplify and then fade to become inaudible.
14. Read at this time the indicated pressure, this reading corresponds to the minimum pressure or diastolic pressure.
15. Unscrew the decompression wheel further to completely purge the remaining air in the bladder.
16. Remove the cuff.

ESSENTIEL  
PERFORMANCE

The sphygmomanometer allows to measure the non invasive blood pressure with a cuff and a stethoscope

## 7. CLEANING, STORAGE AND DISPOSAL

### Cleaning

#### → Manometer cleaning

- A** Spraying of decontaminant solution.  
Do not use detergent, do not immerse

#### → Soft Cuff® & CLINIC® cleaning:

*Recommended product : DETERQUAT / ANIOS.*

- A** Spray the cuff with a detergent, disinfectant and bactericide foam, taking care to spread the product properly.

- B** Leave on for 15 minutes

- C** Rinsing is useless, wipe with a single use tissue

#### → CLINIC® immersion cleaning:

*Recommended product : ANIOS'CLEAN EXCEL D / STERANIOS*

- A** Seal the tube(s) with the closure system provided: • Connector with cap.

- B** Immerse the tube and cuffs in a solution.

- C** Immersion time: 15 minutes.

- D** Rinse thoroughly with water, renewal of the soaking bath after each use.

### Secured disposal of the device

The device should be disposed off, after being decontaminated in compliance with current local practices.

It should be discard in a specific disposal container.

### Usage conditions

| TEMPERATURE      | HYGROMETRY                       | PRESSURE            |
|------------------|----------------------------------|---------------------|
| + 10°C to + 40°C | 15% à 85% (without condensation) | 700 hPa to 1060 hPa |

## 8. INFORMATION IN CASE OF SERIOUS INCIDENT

Any serious incident with the device should be reported and notified to the manufacturer and to the competent Health Authority of the State Member where the patient lives.

## 9. WARRANTY

|                       |                                                                                                                                                                                    |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>WARRANTY</b>       | Spengler replaces or repairs any non-functioning device as notified in the manual. The sphygmomanometer is guaranteed for 3 years.                                                 |
| <b>WARRANTY LIMIT</b> | the product has to be used in compliance with the instructions for use and with the foreseen indications . It should not be modified neither accidentally deteriorated before use. |

## 10. SYMBOLES USED

| SYMBOL                                                                             | INFORMATION        | SYMBOL                                                                             | INFORMATION                  |
|------------------------------------------------------------------------------------|--------------------|------------------------------------------------------------------------------------|------------------------------|
|   | Manufacturer       |   | Cuff                         |
|   | Medical Device     |   | Pouch                        |
|   | Reference          |   | CE                           |
|   | Manometer          |   | Consult instructions for use |
|   | Manufacturing date |   | Swiss representative         |
|  | Batch number       |  | User manual                  |

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**SPENGLER SAS**