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Document No. JDMT-4804-108
Article No.: 001
Version: A
Date d'émission: 2021.04

Warnings

- ⚠ Read instructions thoroughly before using digital thermometer.
- ⚠ Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
- ⚠ Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- ⚠ Do not place thermometer battery near extreme heat as it may explode.
- ⚠ Remove battery from the device when not in operation for a long time.
- ⚠ The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- ⚠ Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- ⚠ Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- ⚠ Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- ⚠ After each use, disinfect the thermometer especially in case the device is used by more than one person.
- ⚠ Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- ⚠ Do not use thermometer orally after being used rectally.
- ⚠ For children who are two years old or younger, please do not use the devices orally.
- ⚠ If the unit has been stored at temperatures over 5°C - 40°C, leave it in 5°C - 40°C ambient temperature for about 15 minutes before using it.
- ⚠ Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- ⚠ PORTABLE RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could result.
- ⚠ It is not intended for use in the oxygen rich environment and presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- ⚠ Do not put the thermometer in direct sunlight or with cotton wool, otherwise the accuracy will be affected.
- ⚠ The equipment should not be cleaned and disinfected while in use.
- ⚠ Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/ or patient is established.

INTENDED USE

The predictive digital thermometers are intended to measure the human body temperature in regular mode orally, rectally and under the arm. And the devices are reusable for clinical or home use on people of all ages, including children under 8 years old with adult supervision.

INTENDED USER/INDICATION FOR USE

Patient or healthy people who want to measure their body temperature, and medical workers who want to measure the patient's body temperature.

Patient target group :Patients or healthy people who need or want to have their temperature taken.

PLEASE READ CAREFULLY BEFORE USING

The predictive digital thermometer provides a quick and highly accurate reading of an individual's body temperature. Predictive-read thermometers are quicker than actual read thermometer. Therefore, users only need about 20 seconds to take the temperature readings. Due to measurement sites of predict thermometer are different, the reading time also may be different, but most people's actual time is usually between 12 and 30 seconds.

This appliance conforms to the following standards:

ISO 80601-2-56 Medical electrical equipment – Part 2-56:Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement,
EN 60601-1-11 Medical electrical equipment – Part 1-11: General requirements for basic safety and essential performance


Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of EN 60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

And the manufacturer is ISO 13485 certified.

CONTENTS

1 thermometer - 1 user's manual- 1 Storage Case

PRODUCT ILLUSTRATION

PRECAUTIONS

→ The performance of the device may be degraded should one or more of the following occur:

- Operation outside the manufacturer's stated temperature and humidity range.
- Storage outside the manufacturer's stated temperature and humidity range.
- Mechanical shock (for example, drop test) or degraded sensor.
- Patient temperature is below ambient temperature.

→ Portable and mobile RF communications can affect the device. The device needs special precautions regarding EMC according to the EMC information provided in the accompany documents.

→ Do not use the devices in the MR environment.

SYMBOL EXPLANATION

| Symbol | Description | Symbol | Description |
|--|--|---|---|
|  | Caution ! |  | Date of manufacture |
|  | Refer to instruction manual/ booklet |  | Storage and Transportation Temperature Limit (-20°C - 55°C) |
|  | Type BF applied part |  | Storage and Transportation Humidity limitation (15% - 95% HR) |
|  | Batch Code |  | Atmospheric pressure limitation |
|  | Authorized representative in the European Community |  | Medical device |
|  | Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products. |  | French disposal recommendations |
|  | CE conformity marking 0123 represent Notified Body identification number |  | MR Unsafe |
|  | Manufacturer |  | The first num.2:Protected against solid foreign objects of Ø 12.5 mm and greater. The second num.7:Protected against the effects of temporary immersion in water. |
|  | Imported by |  | Distributed by |
|  | Catalog reference | | |

SPECIFICATIONS

| | | | | |
|--------------------------------------|---|----------|------------|------------|
| Type | Digital Predictive Thermometer | | | |
| Measure Range | 32,0°C - 43,9°C | | | |
| Accuracy | ±0.1°C during 35.5°C-42.0°C at 18°C-28°C ambient operating range ±0.2°C for other measuring and ambient operating ranges | | | |
| Average measure time | 20 seconds | | | |
| Operating mode | Adjusted mode: Predictive mode | | | |
| Display | Liquid crystal display, 3 1/2 digits | | | |
| Memory | For storing the last measured value | | | |
| Battery | One 1.5 V DC. button battery (size LR41 or SR41, UCC 392) | | | |
| Battery life | Approx. 200 hours of continuous operation or 1 year with 3 measurements per day | | | |
| Dimension | 13.8 cm x 2.2 cm x 1.2 cm (L x W x H) | | | |
| Weight | Approx. 12 grams including battery | | | |
| Expected service life | Three years | | | |
| Ambient operating range: | Temperature : 5°C - 40°C - Hygrometry : 15 % - 95 % HR Atmospheric Pressure : 70 kPa - 106 kPa | | | |
| Storage and transportation condition | Temperature : -20°C - 55°C - Hygrometry : 15 % - 95 % HR Atmospheric Pressure : 70 kPa - 106 kPa | | | |
| Clinical precision | | Oral Use | Rectal Use | Armpit Use |
| | Clinical bias | -0,21°C | -0,20°C | 0,04°C |
| | Agreement limits | 0,52°C | 0,47°C | 0,38°C |
| | Clinical repeatability | 0,08 | 0,09 | 0,08 |
| Ingress Protection Rating: | IP 27 | | | |
| Classification | Type BF  | | | |
| Contraindication | No contraindication | | | |

DIRECTIONS

- Press the On/Off Button next to LCD display. A tone will sound as the screen shows (See Figure 2). After showing the last recorded temperature, the thermometer will enter into the testing mode (See Figure 3).

Figure 2


Figure 3

- Position thermometer in desired location (mouth, rectum, or armpit.)
 - Oral Use:** Place thermometer under tongue as indicated by position «√» shown in Figure 4. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air.
 - Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2) into rectum.
 - Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required.
- The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signaling tone (beep) must be maintained without exception. The temperature will enter the state of maximum update after the buzzer notification. So that in order to achieve better body temperature measurement result, recommend to keep the probe in mouth and rectum about 2 minutes, or in armpit about 5 minutes regardless of the beep sound.

Note: Normally the buzzes are "Bi-Bi-Bi-Bi-"; alarm beeps more rapidly when temperature reaches 37.8°C or higher, and the buzzes are "Bi-Bi-Bi-----Bi-Bi-Bi-----Bi-Bi-Bi"
- To prolong battery life, press the On/Off Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 15 minutes.
- Store the thermometer in its protective case.

Figure 4

TROUBLESHOOTING

| Error message | Problem | Solution |
|--|---|---|
| Lo | Temperature taken is lower than 32.0°C | Turn off, wait one minute and take a new temperature via close contact and sufficient rest. |
| Hi | Temperature taken is higher than 43.9°C | Turn off, wait one minute and take a new temperature via close contact and sufficient rest. |
| Err | The system is not functioning properly. | Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service. |
|  | Dead battery: Battery icon is flashing, can't be measurable | Replace the battery. |

BATTERY REPLACEMENT

- Replace battery when the symbol «» in the lower right corner of the LCD display starts flashing.
- Pull battery cover off as shown in Figure 5.
- Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2) (See Figure 6)
- Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41,UCC392, or equivalent. Be sure battery is installed with "+" polarity facing up. (See Figure 7)
- Slide battery chamber back into place and attach cover.

Note : The following schematic diagram of battery replacement is from the back of the thermometer

Figure 5


Figure 6


Figure 7

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used according to the use instruction, periodic readjustment is not required. However, we recommend checking calibration every two years or whenever clinical accuracy of the thermometer is in question. Set the temperature of the high-precision water bath temperature to 33.00°C, 37.00°C, 42.00°C. Turn on the thermometer and insert into the high-precision water bath and then check the laboratory accuracy of thermometer. Repeat the previous measurement process after it is stable for 10 minutes.

Please send the complete device to the dealers or manufacturer.

The above recommendations do not supersede the legal requirements. The user must always comply with legal requirements for the control of the measurement, functionality, and accuracy of the device which are required by the scope of relevant laws, directives or ordinances where the device is used.

CLEANING AND DISINFECTION

- Immerse the thermometer probe in distilled water for at least 1 minute;
- Using a clean, soft cloth to wipe the thermometer down to remove any residue;
- Repeat step 1 and 2 for three times until no soil is seen with visual inspection after cleaning;
- For thoroughly clean and disinfection, please use method A or B:
 - Method A (High level disinfection):** immerse the thermometer probe in 0.55% OPA(O-Phthaldehyde), such as CIDEX OPA, for at least 12 minutes under temperature at 20°C;
 - Method B (Low level disinfection):** Using a clean soft cloth dipped in 70% medical alcohol, wipe the probe 3 times, at least one minute for each time.
- Repeat step 1 to 3 to remove OPA residuals;
 - Note 1 :** Rectal use is not recommended for home use as OPA will not be readily available outside of a hospital. If rectal measurement is necessary, we strongly recommend high level disinfection.
 - Note 2 :** Please operate according to the manual of OPA for reference.

To prevent damage to the thermometer please note and observe the following:

- Do not use benzene, paint thinner, gasoline or other strong solvents to clean the thermometer.
- Do not attempt to disinfect the sensing probe (tip) of the thermometer by immersing in alcohol, OPA or in hot water (water over 50°C) for long time.
- Do not use ultrasonic washing to clean the thermometer.

LIMITED WARRANTY

The thermometer is guaranteed for one year from the date of purchase. If the thermometer does not function properly due to defective components or poor workmanship, we will repair or replace it free of charge. All components are covered by this warranty excluding the battery. The warranty does not cover damages to your thermometer due to improper handling. To obtain warranty service, an original or copy of the sales receipt from the original retailer is required.

 Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.

 The product is in compliance with the requirements of MDR(EU)2017/745, «0123» is the identification number of notify body.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regards to EMC which must be published in the IFU. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1 : Electromagnetic Compatibility Information - Electromagnetic emissions

| The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | |
|--|--------------------|
| Emissions Test | Compliance |
| CISPR 11 Conducted emissions | Not applicable |
| CISPR11 Radiated emission | Group 1 Class B |
| Harmonic emissions IEC 61000-3-2 | Not applicable |
| Voltage fluctuations/flicker emissions IEC 61000-3-3 | Not applicable |

Table 2 : Electromagnetic Compatibility Information - Electromagnetic immunity

| The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | |
|--|---|---|
| Immunity Test | IEC 60601 Test level | Compliance level |
| Electrostatic discharge (ESD) IEC 61000-4-2 | +/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air | +/- 8 kV contact +/- 2 kV, +/- 4 kV, +/- 8 kV, +/- 15 kV air |
| Electrostatic transient/burst IEC 61000-4-4 | +/- 2 kV for power supply lines 100 kHz repetition frequency +/- 1 kV for input/output lines | N/A |
| Surge IEC 61000-4-5 | +/- 0,5 kV, +/- 1 kV differential mode line-line | N/A |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0% UT (100% dip in UT) for 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT (100% dip in UT) for 1 cycle at 0° 70% UT (30% dip in UT) for 25/30 cycles at 0° 0% UT (100% dip in UT) for 250/300 cycles at 0° | N/A |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 30 A/m, 50/60 Hz | 30 A/m, 50/60 Hz |

NOTE : UT is the A.C. mains voltage prior to application of the test level

Table 3 : Electromagnetic Compatibility Information - Electromagnetic immunity

| The device is intended to be used in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. | | |
|--|--|---|
| Immunity test | IEC 60601 Test level | Compliance level |
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 80 MHz outside ISM bands | N/A |
| Radiated RF IEC 61000-4-3 | 10 V/m 80 MHz to 2,7 GHz | 10 V/m |
| Proximity magnetic fields IEC 61000-4-39 | 30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m Pulse modulation 2,1 kHz 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz | 30 kHz, 8 A/m, CW 134,2 kHz, 65 A/m Pulse modulation 2,1 kHz ; 13,56 kHz, 7,5 A/m, Pulse modulation 50 kHz |

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflections from structures, objects and people.

- a. The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are the following: 6 765 MHz to 6 795 MHz ; 13 553 MHz to 13 567 MHz ; 26 957 MHz to 27 283 MHz ; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz et 80 MHz are the following: 1,8 MHz to 2,0 MHz ; 3,5 MHz to 4,0 MHz ; 5,3 MHz to 5,4 MHz ; 7 MHz to 7,3 MHz ; 10,1 MHz to 10,15 MHz ; 14 MHz to 14,2 MHz ; 18,07 MHz to 18,17 MHz ; 21,0 MHz to 21,4 MHz ; 24,89 MHz to 24,99 MHz ; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- b. The compliance levels in ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that portable/mobile communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in this frequency ranges.
- c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephone and land mobile radio, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary such as re-orienting or relocating the device.
- d. Over the frequency range 150 kHz to 80 MHz, field strength should be less than 3 V/m.

Table 4 : Recommended separation distances between RF wireless communication equipment

The device is intended to be used in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interferences by maintaining a minimum distance between RF wireless communications equipment and the device as recommended below, according to the maximum output power of the communications equipment.

| Frequency MHz | Maximum Power W | Distance | Test level IEC 60601 | Compliance level |
|---------------|-----------------|----------|----------------------|------------------|
| 385 | 1,8 | 0,3 | 27 | 27 |
| 450 | 2 | 0,3 | 28 | 28 |
| 710 | 0,2 | 0,3 | 9 | 9 |
| 745 | | | | |
| 780 | | | | |
| 810 | | | | |
| 870 | 2 | 0,3 | 28 | 28 |
| 930 | 2 | 0,3 | 28 | 28 |
| 1 720 | | | | |
| 1 845 | | | | |
| 1 970 | | | | |
| 2 450 | | | | |
| 5 240 | 0,2 | 0,3 | 9 | 9 |
| 5 500 | | | | |
| 5 782 | | | | |

NOTE 1 : These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflections from structures, objects and people.

WARNINGS !

- ➔ This device should not be used in the vicinity or on the top of other electronic equipment, such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- ➔ The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement part for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- ➔ The operator should maintain basic safety

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Document reference:
NU_TEMPOp_1VB_20241112

Document No. JDMT-4804-108
Article No.: 001
Version: A
Date of Issue: 2021.04

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SPENGLER SAS
JOYTECH Healthcare Co., Ltd.
**Shanghai International Holding Corp. GmbH (Europe)**