Spengler

TEMPO Ear Pro+

DET-1215

CE 0123

Thermomètre auriculaire infrarouge
Ear Infrared thermometer
Infrarot-Ohrthermometer
Termómetro de oído por infrarrojos

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European Authorized Representative:

Shanghai International Holding Corp. GmbH
(Europe) Eiffestrasse 80, 20537 Hamburg, Germany

Swiss Mandatary

MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug Switzerland

The product is in compliance with the requirements
of MDR (EU) 2017/745, "0123" is the identification
number of notify body

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INTRODUCTION

Please read all instructions carefully and thoroughly before using this product.

The TEMPO Ear PRO+ thermometer is specifically designed for safe use in the eardrum. The Infrared Ear Thermometer is a device capable of measuring people's body temperature by detecting the intensity of infrared light emitted from the ear canal of human. It converts the measured heat into a temperature reading displayed on the LCD. When properly used, it will quickly assess your temperature in an accurate manner. Intended use: The infrared ear thermometer is intended for the intermittent measurement of human body temperature. The device can be reused by people of all ages for home use and clinical use. Indications for Use: The infrared ear thermometer is used to measure body temperature. This appliance conforms to the following standards:

ISO 80601 2-56 Medical electrical equipment

- Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement,

IEC 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral

Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment and complies with the requirements of IEC 60601-1-2(EMC), IEC 60601-1(Safety) standards. And the manufacturer is ISO 13485 certified.


WARNING

1. It's not meant to replace a visit to the doctor. Please consult with doctor if you have health concerns.
2. There is no gender or age limitation for using infrared ear thermometer.
3. To avoid inaccurate readings, always attach a new, clean probe cover for each temperature reading.
4. Do not expose the thermometer to temperature extremes (below -25°C/-13°F or over 55°C/ 131°F) nor excessive humidity (>95%RH).
5. This thermometer must only be used with this brand thermometer probe covers.


6. Keep the battery away from children.
7. Remove battery from the device when not in operation for a long time.
8. Do not put the thermometer in direct sunlight or with cotton wool, otherwise the accuracy will be affected.
9. Portable and mobile RF communications can affect the device. The device needs special precaution regarding EMC according to the EMC information provided in the accompany documents.
10. ME equipment should not be cleaned while in use.
11. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
12. The probe of the ME equipment shall not be serviced or maintained while in use with a patient.

PRODUCT DESCRIPTION

- | | |
|---|--------------------------------|
| 1. Probe Cover | 8. LCD |
| 2. Probe | 9. Probe Cover ejection button |
| 3. Test Button | |
| 4.  Button | 10. Probe Cover box |
| 5. Setting Button | 11. Battery cover |
| 6. Timer Button | 12. Reminder light |
| 7. Charging base | |

LCD DISPLAY INTRODUCTION

BASIC FUNCTIONS

Real Time Clock	<p>The real time clock will be recorded with the memory function and help you to recognize each measurement result.</p> <p>Please see the Real time clock setting section to learn how to setup the time in the first use.</p>
Ear Temperature Mode	<p>The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature. Please see the Illustration For Use section to learn how to measure the body temperature.</p>
Probe cover detection and eject function	<p>The thermometer has been designed for use with probe cover. Users can press the eject button to remove the probe cover.</p>
Beep Hint	<p>When the testing is done, it will sound one short beep.</p>
Memory mode	<p>The machine can store up to 30 sets of measurements. Each memory also records the measurement results date/time/ mode. When the data exceeds 30 sets, the coverage starts from the earliest data.</p>

°C / °F Switch	Please see the Real Time Clock Setting section to learn how to change between Celsius and Fahrenheit.
Sound Switch	You can use the Setting button to choose whether to turn on or off the sound.
Timer function	The thermometer includes a 60 second timer that features an audible notification and visual indicator at 0, 15, 30, 45, and 60 seconds. Please see the Illustration For Use section to learn how to measure the body temperature.

EAR THERMOMETER ADVANTAGES

Infrared Ear Thermometer measures core body temperature, which is the temperature of a body's vital organs. (See Figure 1) Ear temperatures accurately reflect core body temperatures since the eardrum is in the cranial cavity and is not affected by the environment. The eardrum shares blood supply with the hypothalamus (part of the brain.) Therefore, body temperature changes are reflected sooner in the ear than they are in other parts of the body.


Figure 1

EAR THERMOMETER ADVANTAGES

Clinical studies have shown that the ear is an excellent site for temperature measurement because temperatures taken in the ear reflect the body's core temperature. Body temperature is regulated by the hypothalamus, which shares the same blood supply as the tympanic membrane. Changes in core body temperature are usually seen sooner at the tympanic membrane than at other sites, such as the rectum, mouth or under the arm.

Advantages of taking temperatures at the ear versus other sites:

- Axillary temperature readings only reflect skin temperature which may not indicate the internal body temperature.
- Rectal temperatures often lag significantly behind internal body temperature changes, especially at times of rapidly changing temperatures.
- Oral temperatures are often influenced by eating, drinking, thermometer placement, breathing through the mouth, or the inability of the person to close their mouth completely.



BODY TEMPERATURE

The temperature of a healthy person is affected by various factors: the person's individual metabolism, their age (body temperature is higher in babies and toddlers and decreases with age. Greater temperature fluctuations occur faster and more often in children, c.g. due to growth spurts), their clothing, the ambient temperature, the time of day (body temperature is lower in the morning and increases throughout the day towards evening), the preceding physical and, to a lesser extent, mental activity.



It is suggested that users develop the habit of daily measurement, and use this product to establish a private normal temperature range for themselves and their families as a reference for measurement.

REAL TIME CLOCK SETTING

When using thermometer for the first time, please set the parameters of the thermometer. With the thermometer off, press and hold Setting BUTTON to enter into setting mode.

1. Set the unit Press  button to select the unit you want. After the unit is set, press Setting Button, the time format figure will appear .

2. Set the time format The device can display the time in either an AM/PM (12-hour) or a 24:00 (24-hour) format. Press and release  button to select the format. With the preferred time format on the display, press Setting Button, the Hour figure is flashing automatically.

3. Set the hour Press and release the  button to advance one hour until the correct hour appears. After the hour is set, press Setting Button, the Minute figure is flashing automatically.

REAL TIME CLOCK SETTING

4. Set the minute

Press and release the  button to advance one minute until the correct minute appears. After the minute is set, press Setting Button, the Year figure is flashing automatically.

5. Set the year

Press and release the  button to advance one year until the correct year appears. After the year is set, press Setting Button, the Month figure will appear.

6. Set the month

Press and release the  button to advance one month until the correct month appears. After the month is set, press Setting Button, the Date figure is flashing automatically.


7. Set the date

Press and release the  button to advance one day until the correct day appears. After the day is set, press Setting Button to exit the setting mode.

INSTRUCTIONS FOR USE

To measure ear temperature:

1. To achieve accurate readings, if a probe cover is needed, make sure a new, clean probe cover is in place before each measurement. Place a new probe cover on the probe firmly; making sure the probe and the inner portion of the probe cover fit snugly(Figure 2). If the probe cover needs to be removed after the measurement is completed, press the probe cover ejection button with appropriate force to remove the probe cover (Figure 3).

2. Press the  button, the display is activated to show all segments. After self-checking Figure 4 appears on the display screen with last memory data (Figure 5), then the thermometer automatically enters the ear temperature measurement mode (Figure 6) and three lines appear in the center of the screen successively with ear symbol  and remain stationary, then you can start a new measurement. You can measure temperature again once the green reminder light and ear symbol reappear.

Remarks:

- a. TEMPO Ear PRO+ must be used with probe cover, if not installed with probe cover, then "CAP" appears on the display with  symbol flashes from left to right (as shown in Figure 7)
- b. To avoid the risk of cross contamination, it is recommended that you always measure the ear temperature with a new probe cover.

Figure 4

Figure 5

Figure 6


Figure 7

3. Insert the probe along the ear canal as deep as possible into the ear as depicted in Figure 8. Then press the Test BUTTON once to start the measurement. A long beep will sound approximately 1 second after the Test BUTTON is activated. This signals the end of the process as results are displayed on the LCD screen.
4. The thermometer is ready for testing again once an ear mark displays on the LCD. A waiting period of 20 seconds between testing is recommended to avoid excessive cooling of the skin.
5. The thermometer will shut off automatically after 60 seconds of inactivity. To prolong battery life, press  button to turn the unit off.


Figure 8

How to turn on or off sound:

You can press Setting Button to turn on or off sound.

After measurement:

1. Power off:
 - Device will automatically shut off if left idle for more than 10 seconds to extend battery life.
 - Press  button to turn it off .
2. Clean the probe after each use to ensure an accurate reading and avoid cross contamination.
(See the section of Care and Cleaning for details.)

BackLight:

In Ear mode:

1. The display will be lighted GREEN for 3 seconds with a happy face 😊 when the unit is ready for measurement and a measurement is completed with a reading less than 37.3°C (99.1°F).
2. The display will be lighted YELLOW for 3 seconds with a happy face 😊 when a measurement is completed with a reading less than 37.8°C (100.0°F).
3. The display will be lighted RED for 3 seconds with a bad face ☹️ when a measurement is completed with a reading equal to or higher than 37.8°C(100.0°F).

Reminder light:

1. Blinking:

When the temperature measurement is in progress.

2. Constant light on:

When the temperature measurement is complete and the probe position does not change during the temperature measurement, the unchanged probe placement position helps to measure the accurate temperature.

Pulse auxiliary measurement function:

The infrared ear thermometer includes a 60 second timer that features an audible notification and visual indicator at 0, 15, 30, 45, and 60 seconds (See Figure 9). Press the Timer button to active timer. A beep is heard at initiation of the timer. The display will show timer counting up in seconds. The timer will beep at the completion of each 15 second interval to provide audible notification. Then the current segment will become solid and the next segment will flash. At 60 seconds three beeps will sound, all quadrants will appear solid ending the timer function. The thermometer will exit the timer mode after timing is over. Users can also press the timer button to close the function at any time. This feature can be used to time pulse, respirations, etc.


Figure 9

TEMPERATURE TAKING HINTS

1. The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.
2. The ear must be free from obstructions or excess earwax buildup to take an accurate reading.
3. External factors may influence ear temperatures such as:
 - lying on one ear or the other
 - their ears covered
 - exposed to very hot or very cold temperatures
 - recently swimming or bathing
 In these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature.
4. For persons wearing hearing aids or ear plugs, remove the device and wait 20 minutes prior to taking a temperature.
5. When using the thermometer on infants under age 1, pull the ear up making sure the sensor faces the eardrum. (See Figure 10)
6. When using the thermometer on individuals over the age of 1, pull the ear back making sure the sensor faces the eardrum. (See Figure 11)
7. Do not touch the probe window. A cotton swab and rubbing alcohol may be used to clean the sensor window.
8. If the thermometer is stored in a significantly different environment than testing location, place it in the testing location for approximately 30 minutes prior to use.
9. It is not intended for use in the oxygen rich environment and presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
10. The patient is an intended operator.

Figure 10

Figure 11

MEMORY MODE

1. When the thermometer has been turned on and followed by Figure 6 or finished testing, press and hold the Button for three seconds. The letter MEM will appear at the bottom of the display. (See Figure 12)
2. The thermometer will automatically memorize the last 30 temperature readings. Each memory also records the measurement date/time/mode icon. Each time the Button is pressed, the screen displays past readings that correspond with a number 1-30. The number 1 reflects the most recent reading, while the number 30 reveals the oldest reading stored in memory. (See Figure 13)
3. In the memory mode, symbol will not change. The user can press the Test Button to take new measurements.


Figure 12


Figure 13

CARE AND CLEANING

1. Probe covers are intended for single use only! Absent, damaged, perforated, or soiled probe covers can lead to inaccurate temperature readings.
2. The probe window must be kept clean, dry, and undamaged at all times to ensure accurate readings. The accuracy of temperature readings can be affected by damage to the probe window, or the presence of dirt and ear wax on the probe window.
3. Fingerprints, earwax, dust and other soiling compounds reduce transparency of the window and result in lower temperature readings. If the thermometer is accidentally placed in the ear without a probe cover installed, the window must be cleaned immediately.
4. The probe window is the most delicate part of the thermometer. To safely clean the window, gently wipe its surface with a cotton swab slightly moistened with isopropyl alcohol and immediately wipe dry with a clean cotton swab. After cleaning, allow at least 5 minutes drying time before taking temperatures.

Z **E**

Note: Do not use any chemical other than isopropyl alcohol to clean the probe window. Other cleaning agents could cause damage to the thermometer.

5. Use a soft, dry cloth to clean the thermometer display and exterior.
6. Do not submerge the thermometer. Excess liquid may cause damage to the thermometer.
7. Store the thermometer and probe cover in a dry location, free from dust and contamination and away from direct sunlight.
8. Periodic cleaning of the device following use to prevent patient cross infection.
-Use a soft cloth slightly moistened with a 75% isopropyl alcohol solution to wipe the thermometer and probe. Do not use abrasive cleaners.
9. Ensure that children do not use the instrument unsupervised; some parts are small enough to be swallowed.
10. Do not modify this equipment without authorization of manufacturer.
11. Put the thermometer back to the original packaging after using.

BATTERY CHARGING

1. Charge the product when "  " appears in the bottom of LCD display. (See figure 14).
2. If the battery runs out, this product needs to be charged for approximately 6 hours.
3. Put the product on the charge base. The end of the product will touch the charging contact (figure 15) .
4. Charge the product through a power adapter as shown in figure 16 and the product starts charging.

When the icon    circulates, it indicates that the product is being charged.

When the icon remains in this state , it indicates that the product is being fully charged.

Note: Battery life diminishes with use.


Figure 14



Figure 15


Figure 16

SPECIFICATIONS

Measuring range	32.0°C-43.0°C (89.6°F-109.4°F)
Measuring site	Ear canal(Ear Mode)
Reference body site	Oral (This thermometer converts the ear temperature to display its "oral equivalent.")
Operation mode	Ear mode (Adjust mode)
Laboratory accuracy	Ear mode: $\pm 0.2^{\circ}\text{C}$ (0.4°F) during 35.5°C - 42.0°C (95.9°F - 107.6°F) at 15°C - 35°C (59.0°F - 95.0°F) operating temperature range $\pm 0.3^{\circ}\text{C}$ (0.5°F) for other measuring and operating temperature range.
Display resolution	0.1°C or 0.1°F
Measure time	Approximately three seconds
Operating temperature range:	10°C - 40°C (50°F - 104°F), 15%-85%RH, non-condensing Atmospheric Pressure : 70kPa-106kPa
Storage and transport temperature range	-25°C - 55°C (-13°F - 131°F), 15%-95%RH, non-condensing Atmospheric Pressure : 70kPa-106kPa
Clinical accuracy	0-1 year old: Clinical bias: 0.12°C (0.22°F); Clinical repeatability: 0.11°C (0.20°F); Limits of agreement: 0.39°C (0.69°F) 1- 5 years old: Clinical bias: 0.09°C (0.15°F); Clinical repeatability: 0.11°C (0.20°F); Limits of agreement: 0.37°C (0.66°F) Over 5 years old: Clinical bias: 0.24°C (0.43°F); Clinical repeatability: 0.09°C (0.17°F); Limits of agreement: 0.39°C (0.69°F)
Shock	withstands drop of 3 feet
Dimension	146.4*37.6*57mm
Weight	Approx.107grams(with batteries)
Battery	DC3.7V (1900mAh)
Battery life	Approx. 2200 readings
Applied part	Probe

Expected service life/ Shelf life	Three years
Ingress protecting rating	IP22
Accessories	Battery,probe cover
Contraindications	None

TROUBLESHOOTING

Error message	Problem	Solution
Er 2	The ambient temperature is not within the range between 10°C and 40°C (50°F-104°F).	Place the thermometer in a room for at least 30 minutes at room temperature between 10°C and 40°C (50°F-104°F).
Er 3	The thermometer is placed incorrectly or unsteady.	Read Illustration for use thoroughly and take a new temperature measurement.
Er 4 <small>In calibration mode</small>	The thermometer showing a rapid ambient temperature change.	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: between 10°C and 40°C (50°F-104°F).
Er 5	The thermometer is not functioning properly.	Please contact the manufacturer or the manufacturer's representative.
H ₁	Temperature taken is higher than 43.0°C (109.4°F).	Read Temperature taking hints thoroughly, then check the integrity of the probe cover and take a new temperature measurement.
Lo	Temperature taken is lower than 32.0°C (89.6°F).	Read Temperature taking hints thoroughly, then make sure the probe cover and lens filter are clean, then take a new temperature measurement.




Error message	Problem	Solution
	The thermometer works properly.	Indicates battery is between 70% and 100% of usable battery capacity.
	The thermometer works properly.	Indicates battery is between 30% and 70% of usable battery capacity.
	It indicates that the power is low, but you can continue to measure.	Indicates battery is between 10% and 30% of usable battery capacity. The thermometer will take a proper measurement but batteries must be charged soon.
	The thermometer could not work due to low battery.	Indicates battery is under 10% of usable battery capacity. Charge the battery.

CALIBRATION

The thermometer is initially calibrated at the time of manufacture. If the thermometer is used in accordance with the instructions for use, no periodic calibration or readjustment is necessary. Should there be any doubt regarding the device's clinical accuracy, the entire thermometer must be returned to the retailer or directly to the manufacturer.

Only the manufacturer is authorized to perform calibration of the thermometer, following a validated protocol using a calibration tool compatible with the device. Any calibration carried out by the manufacturer will be accompanied by a calibration certificate and a verification report confirming the device's compliance with measurement tolerances. Any spot checks for drift conducted by the thermometer user must be performed using a "black body" type calibration tool. If non-compliance is identified during a spot check, this must be verified through a full calibration at the manufacturer's facility before the thermometer can be declared non-compliant.

SYMBOL EXPLANATION

	Manufacturer		The product conforms to the requirements of the EC Directive MDR(EU) 2017/745 on medical devices.
	Date of manufacturer		Atmospheric pressure limitation
	Importer		Storage and Transportation Humidity limitation: 15%-95%RH
	Distributor		Storage and Transportation Temperature Limit: -25°C-55°C(-13°F-131°F)
	Authorized representative in the European Community		Do not resterilize
	Authorized Swiss Representative		Do not use if package is damaged and consult instructions for use
	Reference		Non-sterile
	Model number		Caution
	Medical Device		Refer to instruction manual/booklet
	TYPE BF APPLIED PART		General symbol for recovery/recyclable
	Batch Code		Disposal of this product and used batteries should be carried out in accordance with the national regulations for the disposal of electronic products.
	Unique device identifier	   Triman for recyclable packaging	
	The first num.2: Protected against solid foreign objects of 12.5 mm and greater. The second num.2: Protection against vertically falling water drops when ENCLOSURE tilted up to 15°		
	Expiration date		

SERVICE

The thermometer has a limited two-year warranty. Do not attempt to disassemble or repair the thermometer by yourself. Should service be required during or after the warranty period you must contact the manufacturer. Repackage the thermometer carefully in its original packaging or securely pack to avoid damage during shipping. Include the original sales slip indicating the date of purchase, a note describing the problem, and your return address. Send the thermometer prepaid and insured. The lay operator or lay responsible organization should contact the manufacturer or the manufacturer's representative:

- for assistance, if needed, in setting up, using or maintaining the thermometer; or
- to report unexpected operation or events.

WARRANTY

Thermometer is warranted by manufacture to be free from defects in material and workmanship under normal use and service for a period of two years from the date of delivery to the first user who purchases the instrument. This warranty does not cover batteries, damage to the probe window, or damage to the instrument caused by misuse, negligence or accident, and extends to only to the first purchaser of the product. Additionally this warranty becomes void, if the thermometer is operated with anything other than this brand thermometer probe covers.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1

Guidance and declaration of manufacturer electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below.		
Emissions test	Compliance	Electromagnetic environment
Radiated emission CISPR 11	Group 1, class B	The device uses RF energy only for its internal function Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	Group 1, class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ±2 kV,±4 kV, ±8 kV, ±15 kV air	± 8 kV contact ±2 kV,±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV , 100kHz, for AC power port	± 2 kV . 100kHz, for AC power port	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5kV, ±1kV (differential mode)	±0.5kV, ±1kV (differential mode)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

<p>Radiated RF EM</p>	<p>3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz</p>	<p>3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz $800 \text{ MHz} < f < 2.7 \text{ GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>
<p>Conducted disturbances Induced by RF fields IEC 61000-4-6</p>	<p>3 V in 0.15 MHz- 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz</p>	<p>3 V in 0.15 MHz- 80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz $800 \text{ MHz} < f < 2.7 \text{ GHz}$ where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: </p>

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity			
<p>Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-type Fully Automatic Digital Blood Pressure Monitor has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.</p>			
Test frequency (MHz)	Band (MHz)	Service	Modulation
385	380-390	TETRA 400	Pulse modulation 18Hz
450	430-470	± 2 kV , 100kHz, for AC power port	GMRS 460 FRS 460
710	704-787	LTE Band 13, 17	Pulse modulation 217Hz
745			
780			
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18Hz
870			
930			
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz
1945			
1970			
2450	2400-2570	Bluetooth,WLAN, 802.11 b/g/n,RFID 2450,LTE Band 7	Pulse modulation 217Hz
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217Hz
5500			
5785			

Table 3

Guidance and declaration of manufacturer-electromagnetic immunity		
<p>Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. Arm-type Fully Automatic Digital Blood Pressure Monitor has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipment and this medical equipment and/or systems as recommended below.</p>		
Maximum power (W)	Distance (m)	Immunity test level (V/m)
1.8	0.3	27
2	0.3	28
0.2	0.3	9
2	0.3	28
2	1700-1990	28
2	2400-2570	28
0.2	0.3	9

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device.		
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Table 5

Test specifications for ENCLOSURE PORT IMMUNITY to proximity magnetic fields.		
Test Frequency	Modulation	IMMUNITY TEST LEVELS (A/m)
30 kHz a)	CW	8
134,2kHz	Pulse Modulation b) 2,1kHz	65 c)
13,56MHz	Pulse Modulation b) 50kHz	7,5 c)
<p>a) This test is applicable only to ME EQUIPMENT and ME SYSTEMS intended for use in the HOME HEALTHCARE ENVIRONMENT.</p> <p>b) The carrier shall be modulated using a 50% duty cycle square wave signal.</p> <p>c) r.m.s., before modulation is applied.</p>		

ELECTROMAGNETIC COMPATIBILITY INFORMATION

WARNINGS !

- This device should not be used in the vicinity or on the top of other electronic equipment such as cell phone, transceiver or radio control products. If you have to do so, the device should be observed to verify normal operation.
- The use of accessories and power cord other than those specified, with the exception of cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment or system.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Portable and mobile RF communications can affect the device.

The device needs special pre-cautions regarding EMC according to the EMC information provided in the accompany documents.

- Do not use the devices in the MR environment.
- The Operator should not use the system and should inform the customer service, if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES.
- PRECAUTION: The performance of the device may be degraded should one or more of the following occur:
 - Operation outside the manufacturer's stated temperature and humidity range.
 - Storage outside the manufacturer's stated temperature and humidity range.
 - Mechanical shock (for example, drop test) or degraded sensor.
 - Patient temperature is below ambient temperature.

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JOYTECH Healthcare Co., Ltd.