



Spengler

SINCE 1907



Magister II

SINGLE CHESTPIECE
STETHOSCOPE FOR A VERSATILE AUSCULTATION



Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr

Tel : +33 (0)4 42 52 89 80
Mail : export@spengler.fr

STETHOSCOPE FOR A VERSATILE AUSCULTATION

Magister II SINGLE CHESTPIECE

Professional.

Better results, more versatile, more professional.



STAINLESS STEEL CHESTPIECE

Stainless steel offers the best acoustic properties, so it was selected for use in the chestpieces for the MAGISTER II stethoscopes. In addition, we increased the thickness and the polishing of this chestpiece in order to refine the acoustics.



ADAPTED TO CUFFS

The MAGISTER II stethoscopes are designed with thin edges, perfectly slipping under the tensiometer cuff.



FLOATING MEMBRANE

Present on both chestpiece:

- Light pressure to hear low frequencies
- Firm pressure to hear high frequencies



DESIGNATION	REFERENCE
● Magister® II stethoscope Adult single chestpiece Silver	507 250
● Magister® II stethoscope Adult single chestpiece Black	507 260

SPENGLER QUALITY

CE mark

Spengler stethoscopes, including Magister II are class I medical devices. They are compliant with European Regulation 2017/745 and with guidelines for safety, health, hygiene and protection of the environment.

After sales service

Spengler SAS provides maintenance and repair of equipment. Go to the Spengler website in the Contact section.

Spengler Warranty

Spengler guarantees the stethoscope Magister II for 5 years from date of purchase. The guarantee does not apply to damage caused by improper handling,

accidents or alterations made to the instrument by third parties. Spengler guarantees the ear tips, diaphragm and non-chill ring for 1 year from date of purchase.

The guarantee does not apply to fair wear and tear of the product/damages and breakdowns due to shock, misuse or neglect/any product disassembled, repaired or modified by the user or unauthorised persons.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

Tel : +33 (0)4 42 52 89 80

Mail : export@spengler.fr