

75



Drap de survie/Rescue Blanket

ISOTHERMIQUE - ISOTHERMAL

Film polyester métallisé très résistant et imperméable.
Réflexion à 90% du rayonnement infrarouge.

Very strong and waterproof metallized polyester film.
Reflects 90% of the infrared radiations.

INSTRUCTIONS D'UTILISATION ET MISE EN GARDE :


La face dorée placée à l'extérieur permet de conserver la température du corps et garantit une meilleure protection contre l'humidité, le froid et les intempéries. La face argentée placée à l'extérieur protège de la chaleur et de l'insolation. Ne pas utiliser pendant un orage, lors d'une intervention chirurgicale si un appareil d'électrochirurgie est utilisé, à proximité de matières incandescentes ou sur la plage.


100

INSTRUCTIONS OF USE & CAUTION:

The golden side on the outside maintains the body temperature and guarantees a better protection against moisture, cold and bad weather. The silver side placed on the outside protects from heat and sunstroke. Do not use during storm, medical surgery if electrical equipment is used and close to incandescent substance or on the beach.

2 faces/sides: **argent/silver - or/golden**
Dimension/size: **220 x 140 cm.**
Poids/Weight: **57 g.**
Non-stérile/non-sterile.

 Importé et distribué par Holtex
30 rue Jean de Guiramaud- 13290
Aix-en-Provence - France

 Wellman Healthcare Tech. (Suzhou) Co., Ltd.
1-3603, Global 188, No.88 Suzhou Road, Suzhou Industrial Park, Jiangsu,
China, 215021

 CMC Medical Devices & Drugs S.L.
C/ Horacio Lengo Nº 18, CP 29006, Málaga, Spain



MD REF SCOUV03

