

# DISH. KERN WEIGHING MWA 300K-1PM FOLDING BAR. DOOR POSTER. CLASS III

KERN

Ref. EKMWA\_1PM\_CIII





CE 0297 1/2



30/08/2024 08:26:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France www.gsh-med.fr **↓** +33(0)4 42 90 31 31**⋈** info@gsh-med.fr



# DISH. KERN WEIGHING MWA 300K-1PM FOLDING BAR. DOOR POSTER. CLASS III



2/2

#### Ref. EKMWA\_1PM\_CIII

### Wheelchair weighing platform with folding barrier and display door approved

## **Description**

- · With folding barrier and display door
- . The large steel weighing pan is also ideal for weighing hard, overweight patients
- . Stable and non-slip thanks to its height-adjustable rubber feet
- · Hold function: for agitated patients, the average value looks for a stable weighing value and "freezes"
- · BMI function for calculating underweight/normal/overweight
- · Large LCD display with backlight, height of figures 25 mm
- · Graduation 0.1 kg
- · Without wheels
- · Battery operated (6 x 1.5V AA not included) or mains (power supply supplied) or optional with rechargeable battery pack
- · Approved in class III
- · Medical approval according to 93/42/EEC

Technical specifications:

· Display dimensions :L×

P×H 200×128×55 mm

- $\cdot$  Total dimensions, without restraint bar: L×P×
- H 1180×1000×1130 mm
- · Display cable length approx. 180 cm
- · Net weight approx. 76 kg

For any measuring device requiring medical Class III approval, accessories must be ordered at the same time as the device so that the approval can be valid.

## **Features**

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040201061	HS Code	84231090
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	90
Bare Unit / Weight (kg)	76		
SPECIFIC STORAGE CONDI	TIONS		
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 08:26:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France www.gsh-med.fr +33(0)4 42 90 31 31 info@gsh-med.fr