

WHEELCHAIR/ SCALE MWA(300KG) WITH FOLDING SIDEBAR/DISPLAY **HOLDER**



Ref. EKMWA_1PM















30/08/2024 08:24:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France







info@gsh-med.fr



www.gsh-med.fr



WHEELCHAIR/ SCALE MWA(300KG) WITH FOLDING SIDEBAR/DISPLAY HOLDER



Ref. EKMWA 1PM

Wheelchair weighing platform with folding barrier and display door

Description

- · With folding barrier and display door
- The large steel weighing pan is also ideal for weighing hard, overweight patients
- . Stable and non-slip thanks to its height-adjustable rubber feet
- Hold function: for agitated patients, the average value looks for a stable weighing value and "freezes"
- · BMI function for calculating underweight/normal/overweight
- · Large LCD display with backlight, height of figures 25 mm
- · Graduation 0.1 kg
- · Without wheels
- · Battery operated (6 x 1.5V AA not included) or mains (power supply supplied) or optional with rechargeable battery pack
- Homologable in Class III
- · Medical approval according to 93/42/EEC

Technical specifications:

· Display dimensions :L×

P×H 200×128×55 mm

- · Total dimensions, without restraint bar: L×P×
- H 1180×1000×1130 mm
- Display cable length approx. 180 cm
- Net weight approx. 76 kg

For any measuring device requiring medical Class III approval, accessories must be ordered at the same time as the device so that the approval can be valid.

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	3661040201054	HS Code	84231090
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	100
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	113
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	118
Packed unit / Weight	90	Bare Unit / Weight (kg)	76
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
acked unit / Quantity 39 - For information purposes only, can be subject to change			

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr



+33(0)4 42 90 31 31



