

EASY FUSE® PRESSURE INFUSION CUFFS

TECHNICAL SHEET

All our pressure infusion cuffs are manufactured in France, in our Manufacture in Castelculier (near Agen)!



SINGLE-USE

With their simple, robust design, these single-use cuffs are economical, while maintaining their safety and reliability:

- High-quality technical fabric.
- Pressure indicator to monitor pressure.
- Transparent veil to visualize the liquid.





Letter	Description	Feature
A	Cuff alone 500 cc	PU-coated fabric + polyamide monofilament
A	Cuff alone 1000 cc	PU-coated fabric + polyamide monofilament
A	Cuff alone 3000 cc	PU-coated fabric + polyamide monofilament
B	Pressure indicator	ABS
C	Clamp	Polyoxymethylene
D	Tubing	PVC
E	3-way tap	Polycarbonate + Polyethylens
F	Bulb	PVC phtalates free

CLEANING & DISINFECTION	Single-use device	
STORAGE CONDITIONS	Packaging type	Origin
	Local	Ventilated
	Temperature	+ 10° C to + 70° C
	Hygrometry	10 to 100 %
	Atmospheric pressure	500 to 1060 hPa

Description	Size	Bladder dimensions (mm)		Reference
		Lenght	Width	
Single-use pressure infusion cuff	500 cc	250	162	M10500
	1 000 cc	365	170	M11000
	3000 cc	338	235	M13000

FRENCH QUALITY DDM

CE mark

Pressure infusion cuffs are Class I medical devices. They comply with current European regulations.

Customer Technical Support Centre

GSH provides pressure infusion cuffs maintenance and repair services. Contact our After-Sales Service on: +33 4 42 90 31 31.

Warranty

DDM pressure infusion cuffs are guaranteed for 24 months (under normal conditions of use). The warranty does not cover faulty maintenance by the user.