

# PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM



### Ref. GMBU02









27/03/2025 21:11:01 - For information purposes only, can be subject to change





## PNEUMATIC MANUAL TOURNIQUET - SINGLE USE - DDM



#### Ref. GMBU02

Dedicated to emergency situations, the single-use DDM manual pneumatic tourniquet is completely stand-alone, making it a true French innovation!

### **Description**

Designed to accompany the patient from the accident site to the emergency room, DDM's single-use pneumatic tourniquet is accessible and easy to use, thanks to its safety valve.

This single-use pneumatic tourniquet is self-sufficient in its composition: bulb, cuff and valve. The injured patient leaves and keeps the tourniquet on until surgery.

With a wide cuff and foam reinforcement inside, pressure is distributed more widely, avoiding skin lesions and intense pain.

Without a manometer, thanks to a mechanical valve, this manual pneumatic tourniquet avoids overpressure, so pressure does not exceed 300mmHg.

This tourniquet also features a label directly on the cuff for recording patient data, such as patient name, date and time of tourniquet application.

#### - Available in :

Adult 23-33cm (GMBU04) Child (GMBU02) XL Adult 38-50cm (GMCU05)

#### **Features**

CE Marking	CE	Medical device	Class I
Bar code	3700468209287	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS	)		
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

27/03/2025 21:11:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand 13290 Aix-en-Provence | France









www.gsh-med.fr