



Spengler

NEW

TEMPO Ear Pro+

The ear infrared thermometer for healthcare professionals in private practice, clinics and hospitals, with Spengler quality.



USAGE HOSPITALIER
CLINIC USE



Spengler

SPENGLER MEDICAL | 30 rue Jean de Guiramaud | 13290 Aix en Provence | FRANCE

CE 0123

www.spengler.fr

TEMPO Ear Pro+

INFRARED EAR THERMOMETER



The Tempo EAR PRO + by Spengler is an infrared ear thermometer. The cover-probe charging and ejection system ensures hygienic and safe use, reducing the risk of cross-contamination.

The temperature measurement is performed accurately and promptly, guaranteeing optimal comfort for both the patient and the practitioner.

It is an essential ally of hospitals, clinics, and medical offices due to its ergonomic, intuitive, and designed for intensive use in the medical environment. The efficiency on a daily basis for reliable and simplified temperature recording.



PRE-HEATING THE PROBE
TO 36°C



STORAGE AND RECHARGE
BASE



PROBE COVERS
CHARGING SYSTEM



MEASUREMENT TIME:
LESS THAN 3 S

30
MEM

MEMORY OF THE LAST
30 TIME-STAMPED
MEASUREMENTS



COUNTER OF
PULSE OF 60S



STANDBY MODE
AFTER 60 S OF INACTIVITY

2 200

BATTERY LIFE OF 2,200
MEASURES



LANYARD INCLUDED



- Measurement time: less than 3 seconds
- Despite movement, the probe's stability ensures accurate temperature measurements.
- Ear temperature range: 32.0 °C to 43 °C (90°F to 109.4 °F)
- Power supply charger included

ALSO AVAILABLE

TEMPO Ear Pro

- Probe cover charging system
- Measurement time: less than 3 s
- Storage base
- **Battery power**



THE SPENGLER QUALITY

CE marking

European Directive EEC 93/42: infrared ear thermometer class IIa. Spengler's ear infrared thermometers are compliant with the vital demands of safety, health, hygiene, and environmental protection.

Spengler Guarantee

Spengler 2-year warranty for TEMPO Ear Pro + ear thermometers from Spengler (under normal conditions of use): free replacement of recognized defective parts by standard exchange (functional defect resulting from a material, manufacturing or design defect) by our after-sales service. Wear parts (batteries are not covered by the warranty). The warranty does not cover maintenance defects on the part of the user.

After-sales service

SPENGLER MEDICAL ensures the maintenance and repair of equipment.

DESIGNATION

REFERENCE

Tempo Ear Pro + Thermometer

THE003

Tempo Ear Pro Thermometer

THE004

Box of 800 probe covers for Tempo Ear Pro Thermometer

THE003ACC1

Rechargeable battery for Tempo Ear Pro + Thermometer

THE003ACC3

Charging station for Tempo Ear Pro Thermometer

THE003ACC2



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr