

NEW



Spengler

TEMPO Ear Pro

The auricular infrared thermometer for healthcare professionals in private practice, clinics and hospitals, with Spengler quality.



CLINIC USE

TEMPO Ear Pro

Spengler

CE 0123

 Spengler

SPENLER MEDICAL | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr

TEMPO EarPro

INFRARED EAR THERMOMETER



The Tempo EAR PRO by Spengler is an ear thermometer with infrared technology. Thanks to its charging and ejection system, it guarantees a hygienic and safe use reducing the risk of cross contamination.

The temperature measurement is accurate and fast ensuring optimal comfort for both patient and practitioner.

Ergonomic, intuitive and designed for intensive use in the medical environment, it is an essential ally of hospitals, clinics and medical practices. Everyday efficiency for a reliable and simplified temperature measurement.



PROBE PREHEATING
TO 36°C



STORAGE BASE



PROBE COVER
CHARGING SYSTEM



MEASUREMENT TIME:
LESS THAN 3 S

30
MEM

MEMORY OF THE LAST
30 MEASUREMENTS



READY TO USE: 2 X AA
BATTERIES AND A REFILL
OF 20 PROBE COVERS
INCLUDED



AUTOMATICALLY
SHUTTING DOWN AFTER
60 S OF INACTIVITY

ALSO AVAILABLE

TEMPO EarPro+

- 60 S pulse counter
- Tie lanyard
- Rechargeable battery
- Charging station



- Measurement time: less than 3 seconds
- Sensor stability feature during temperature measurement for reliable measurement despite movement
- Ear temperature range: 32.0°C to 43°C (90°F to 109.4°F)
- Battery power: 2 AA

MODEL

REFERENCE

Thermometer Tempo Ear Pro	THE004
Thermometer Tempo Ear Pro +	THE003
Box of 800 probe covers for Thermometer Tempo Ear Pro	THE003ACC1
Rechargeable battery for Thermometer Tempo Ear Pro +	THE003ACC3
Charging station for Thermometer Tempo Ear Pro	THE003ACC2

SPENGLER QUALITY

CE marking

European Directive EEC 93/42: auricular infrared thermometer class IIa. Spengler's ear infrared thermometers meet the essential requirements of safety, health, hygiene and environmental protection.

Spengler Guarantee

Spengler 2-year warranty on Spengler TEMPO Ear Pro ear thermometers (under normal use): free replacement of defective parts through standard exchange (malfunction arising from a defect in material, manufacture or design) by our after-sales service. Wear parts (batteries are not covered by warranty). The warranty does not cover maintenance defects on the part of the user.

After-sales service

SPENGLER MEDICAL provides maintenance and repair of equipment.



www.spengler.fr

SPENGLER MEDICAL | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr