



# Spengler

## *Lian® Dessillons*

SPHYGMOMANOMETER

**Spengler continues his story:  
the French production!**



By manufacturing the Lian® Dessillons at its MANUFACTURE in Castelculier (Agen),  
**SPENGLER MEDICAL** definitively anchors its history in French reindustrialization.



CE  
0459



Spengler

SPENGLER MEDICAL | 30 rue Jean de Guirmand | 13290 AIX-EN-PROVENCE | FRANCE

[www.spengler.fr](http://www.spengler.fr)

Tel. : +33 (0)4 42 90 31 31  
Mail : [contact@spengler.fr](mailto:contact@spengler.fr)

# SPHYGMOMANOMETER

## Lian® Dessillons



### A NEAT DESIGN: ERGONOMIC AND ELEGANT

The body design has been rounded to improve the grip.  
The decompression valve is ergonomic and aesthetic.  
The display size allows an easy reading.



### ROBUST AND PRECISE

Its aluminium body, its ABS housing, its TPU protector and its shatterproof glass made it a shock-resistant device.  
Precision clock making mechanism (german made) provides a graduated precision from 0 to 300mmHg with  $\pm 3$ mmHg on the measurement range.



### NEAT QUALITY

Manufactured in France. German clock making mechanism. German bulb.



### FITTED CUFFS

The Soft Cuff® with bladder are comfortable and soft.  
The Clinic® cuffs are adapted to hospitals and clinics as they are immersible for decontamination.



#### SPHYGMOMANOMETERS

DESIGNATION	SOFT CUFF®	CLINIC® CUFF
Spengler Sphygmomanometer Lian® Dessillons, Small adult	TS50003S	TI50003S
Spengler Sphygmomanometer Lian® Dessillons, Adult	TS50004S	TI50004S
Spengler Sphygmomanometer Lian® Dessillons, Adult large	TS50005S	TI50005S
Spengler Sphygmomanometer Lian® Dessillons, Multi cuffs	TS53000S	TI53000S

#### CUFFS ONLY

DESIGNATION	SOFT CUFF®	CLINIC® CUFF
Newborn size bayonet cuff	BT101S	BM101S
Child size bayonet cuff	BT102S	BM102S
Little adult size bayonet cuff	BT103S	BM103S
Adult size bayonet cuff	BT104S	BM104S
Large adult size bayonet cuff	BT105S	BM105S
Lower limb adult size bayonet cuff	CT104S	CM104S

## SPENGLER QUALITY

### CE certificate

Spengler sphygmomanometers as Lian® Dessillons are medical class Im. They comply with European current regulations.

### Customer technical support center

Spengler SAS ensures the maintenance and repair of equipment. Go to the Spengler website in the customer service section or contact the: +33(0)4 42 90 31 31

### Spengler warranty

Spengler 3-year warranty (under normal conditions of use): repair or replacement free of charge of defective parts recognized (malfunction resulting from a defect in material, manufacturing or design) by our after-sales service. The armbands are guaranteed for 12 months. The warranty does not cover maintenance defects on the part of the user.



[www.spengler.fr](http://www.spengler.fr)

SPENGLER MEDICAL | 30 rue Jean de Guiramand | 13290 AIX-EN-PROVENCE | FRANCE

Tel. : +33 (0)4 42 90 31 31

Mail : [contact@spengler.fr](mailto:contact@spengler.fr)