

Distributed by:



NC6 & NC7 COMEN

The NC6 and NC7 multi-parameter monitors: excellence in monitoring!



NC6 & NC7 COMEN

The NC6 and NC7 multi-parameter monitors: excellence in monitoring!

Distributed by:
SPENGLER
MEDICAL



The NC6 and NC7 monitors are adaptable and can be employed in any clinical setting. They are indeed ideal during the service visit, at emergency triage as well as in simple monitoring at the patient's bed. Whether it's routine check-ups or intensive care, they meet all needs in all clinical settings.

ACCURATE AND EFFECTIVE, THEY MEASURE THE FOLLOWING VITAL SIGNS:



NIBP

- Quick measurement: 15s
- Average of measurements
- Measurement during inflation (iNIBP)
- Measurement of orthostatic tension

TEMPERATURE



- Ear thermometer TD-1261
(other possible optional choices)

SpO₂



- All types of patients: infants included
SpO₂ Nellcor and Masimo as an option
- Measurement of pulse rate (FP)
- Measurement of the perfusion index (IP) value
of pulsatile blood flow
- Measurement of breathing.



RESPIRATORY FUNCTION - VENTILATION

- CAPNOGRAPHY: And CO₂ (NC7 only)
- Comen Sidestream
- Masimo Sidestream

ERGONOMIC AND EASY TO USE:

- Touch screen
- Customizable interface
- different possible modes: spot check or monitoring
- alarms
- Barcode scanner

The NC6 and NC7 provide a **advanced clinical diagnostic assistant function** that enhances patient management and healthcare professionals' work :

- The patient's physiological information is entered manually, including their height, weight, blood sugar, O₂ source, and O₂ supplementation, etc...
- Early Alert Score Calculation: MEWS, NEWS2, and PEWS.
- Glasgow coma scale (GCS)
- Screening for Critical Congenital Heart Disease (CCC)
- Pain scoring system with different assessment methods.
- Possibility to set a target objective on a parameter.



SEAMLESS INTEGRATION:

Workflow simplification for healthcare professionals. The NC6 and NC7 integrate perfectly with (CMS) Central Monitoring as well as with SIH or DEP using the HL7 protocol.

PACK CONTENT :

- Main Unit
- An adult NIBP cuff 25-35cm
- An adult SpO2 sensor
- A NIBP tube (2m)
- A power supply cable
- An instruction manual
- A thermometer (optional)

NC6 and NC7 are COMEN monitors proposed by **SPENGLER MEDICAL**
Medical devices of class IIb CE0123 compliant with the European regulation in force relating to medical devices.

Mandatory service maintenance every 2 years.

Manufactured by Shenzhen Comen Medical Instruments Co, Ltd
Distributed by **SPENGLER MEDICAL** - 30 rue Jean de Guiramaud - 13290 Aix-en-Provence - France
www.sepengler-med.fr/ Contact: +33(0)4 42 90 31 31 - info@spengler-med.fr