

BLOOD PRESSURE MONITOR AUTOTENSIO® WRIST

SINCE 1907

The Autotensio® Wrist is the perfect Spengler blood pressure monitor for immediate, easy and quick self-diagnosis. Everywhere and for everyone, it offers reliability and accuracy at home or on the move.

High measurement accuracy

- WHO indicator

Easy to use

- Fully automatic
- 2 buttons
- The bracelet covers a wrist circumference of 13.5 to 21 cm)
- 120 stamped memories
- With a very convenient storage box.
- Comes with 2 AAA batteries



COMPATIBLE WITH:

suivi HTA

The application from the FRHTA
(French Research Foundation for Hypertension)



CE
0459



Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr

Tel. : +33 (0)4 42 52 89 80

Mail : export@spengler.fr

BLOOD PRESSURE MONITOR

AUTOTENSIO® WRIST

TECHNICAL CHARACTERISTICS

Measurement method	Oscillometric	
Response time	less than 1 minute	
Cuff size	From 13,5 to 21,5 cm	
Screen	Display	LCD screen (45 x 30 mm)
		Date and time
		Measurements
		Memories
		Low batteries indicator
Arterial pressure	Measuring range	From 0 to 300 mmHg
	Resolution	1 mmHg
	Precision	+/- 3 mmHg
Pulse	Measuring range	From 30 to 180 beats per minute (BPM)
	Resolution	1 BPM
	Precision	+/- 5 % of the reading
Memory	120 time stamped memories	
Functions	WHO indicator	
	Low batteries indicator	
	Automatic power-off	
Buttons	Start-stop and memory	
Air circuit	Inflation	Micro pump
	Deflation	Solenoid valve
Power supply	Type	2 AAA batteries or Adapter in option
	Autonomy	About 200 measurements
Automatic shutdown	Yes	
Conditions of use	Operating temperature	From 10°C to ~40°C
	Storage condition temperature	From -25°C to ~70°C
	Relative humidity	≤ 93% (without condensation)
Size and dimensions	Dimensions	79 x 66 x 28 mm ((LxWxD)
	Weight	228 g (without batteries)
SPENGLER WARRANTY	2 years (in normal use conditions)	

DESCRIPTION	REFERENCE
● Blood pressure monitor Autotensio® wrist	SPG340-
● Blood pressure monitor Autotensio® wrist - Coral	SPG345
● Display of 6 Autotensio® wrist blood pressure monitors	SPG340P
● Display of 6 Autotensio® wrist blood pressure monitors - Coral	SPG345P

SPENGLER QUALITY

CE certificate

Spengler BPM as Autotensio are class IIa medical devices. They comply with European current regulations.

Customer technical support center

Spengler SAS provides devices maintenance and repair services.
Contact our Customer technical support center on: +33 4 42 90 31 31

Spengler warranty

Autotensio BPM is guaranteed for 2 years (under normal conditions of use). Wearing parts (batteries) are not covered by the warranty. The warranty does not cover faulty maintenance by the user.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

Tel. : **+33 (0)4 42 52 89 80**

Mail : **export@spengler.fr**