



TEMPOEASY

NON-CONTACT INFRARED THERMOMETER

The Tempo Easy takes temperature using infrared and contactless rays.

Very high precision measurements

- Infrared technology
- Instant reading

Easy to use & practical

- Non-contact (distance 3-5 cm)
- 32 memories
- 3 reading modes: body, surface and ambient

Additional Information

- Supplied with two LR6 batteries (AA)

THE ACCURATE, MULTIFUNCTION THERMOMETER

Thanks to its high reliability and 3-function mode, SPENGLER's TEMPO EASY thermometer is ideal in many cases.

Family use: Checking of the temperature of the baby bath.
Using its «surface» mode, take water temperature in less than one second.

Health institution: Checking of ambient temperature.
In EPHAD type institutions, checking ambient temperature is necessary information for the well-being of patients. Using its «ambient» function, you can instantly check the room temperature of each patient's room.

Usage by professionals: Checking ambient temperature.
Thanks to its instantaneous temperature measurement and its light weight, the TEMPO EASY thermometer is the ideal tool for liberal professions.
In addition, the infrared technology means probe covers are not required.



CE
0123
Classe IIa



Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

www.spengler.fr

Tel. : +33 (0)4 42 52 89 89

Mail : contact@spengler.fr

TEMPOEASY

NON-CONTACT INFRARED THERMOMETER

TECHNICAL CHARACTERISTICS

Type	Non-contact thermometer	
Measurement	Mesure infrarouge	
Screen	Display	Backlit LCD screen
		Temperature modes (body, surface and ambient)
		Battery charge icon
		Readings
		Memory
		Sound signal
Temperature	Measuring range	Body : 32.0°C ~ 43.0°C Surface : 0°C ~ 60°C
	Resolution	0.1°C
	Accuracy	32.0°C ~ 34.9°C : ± 0.3°C
		35.0°C ~ 42.0°C : ± 0.2°C 42.1°C ~ 43.0°C : ± 0.3°C
Measuring distance	3 cm ~ 5 cm	
Memory	32 sets	
Unit of measurement	Celsius ou fahrenheit	
Alarme	If temperature is too high (>38°C)	
Batteries	DC 3V (2 pcs AA batteries)	
Consumption	≤ 300 mW	
Autonomy	~ 40 000 mesures	
Automatique power off	After 30 secondes	
Normal using conditions	Ambient temperature	10°C ~ 40°C
	Relative Humidity	≤ 85%
Dimensions and size	Unit size	175 x 50 x 52 mm (LxIxH)
	Unit weight	148 g (including batteries)
SPENGLER WARRANTY	2 years (under normal use conditions)	

DESIGNATION	REFERENCE
Tempo Easy bleu	SP3065
Tempo Easy grey	SP3066
Présentoir de 6 Tempo Easy bleus	SP3065P
Présentoir de 6 Tempo Easy grey	SP3066P

SPENGLER QUALITY

CE Marking

Non-contact thermometer TEMPO EASY Spengler is a class IIa medical device. It complies with current regulations by perfectly meeting the general safety and performance requirements.

Customer technical support center

SPENGLER MEDICAL provides devices maintenance and repair services.
Contact our customer technical support center on +33 4 42 90 31 31 or click on SAV tab on www.Spengler.fr

Spengler warranty

Spengler 2-year warranty (under normal conditions of use): repair or replacement free of charge of defective parts recognized (malfunction resulting from a defect in material, manufacturing or design) by our after-sales service. The armbands are guaranteed for 12 months. The warranty does not cover maintenance defects on the part of the user.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

Tel. : **+33 (0)4 42 52 89 89**

Mail : **contact@spengler.fr**