

NON-CONTACT INFRARED THERMOMETER

TEMPO PRO



Forehead measure



High measuring accuracy



Quick measure



Easy to use



CE
0123
Class IIIa



Spengler

Spengler S.A.S | 30 rue Jean de Guirmand | 13290 Aix en Provence | FRANCE

www.spengler.fr

Tel. : +33 (0)4 42 52 89 80

Mail : export@spengler.fr

NON-CONTACT INFRARED THERMOMETER

TEMPO PRO

TECHNICAL CHARACTERISTICS		
Type	Forehead thermometer	
Measurement	Infrared	
Screen	Display	Backlit LCD screen
	Resolution	0,1°C / F°
Temperature	Measuring range	32,0°C - 43,0°C / 89,6°F - 109,4°F
	Unit of measurement	Celsius or Fahrenheit
	Forehead accuracy	≤ 0,2°C (0,4°F)
Indicators	Low battery indicator	
	Colors indicator	< to 32°C (89,6°F) = yellow Between 32°C and 37,6°C = blue Between 37,6 and 38,4 = yellow > 38,4°C = red
Memories	30 memories	
Batteries	2 x 1,5V AAA batteries	
Normal using conditions	Ambient temperature	16,0°C - 35,0°C (60,8°F - 95,0°F)
	Relative humidity	≤ 85%
	Atmospheric pressure	700 hPa - 1060 hPa
Normal storage conditions	Ambient temperature	-25,0°C - + 55,0°C (-13,0°F - +131°F)
	Relative humidity	≤ 95%
	Atmospheric pressure	500 hPa - 1060 hPa
Dimensions and weight	Unit size	113 x 51 x 162 mm
	Unit weight	130 g (avec les piles)
SPENGLER WARRANTY	2 years (under normal use conditions)	

DESIGNATION	REFERENCE
Forehead thermometer TEMPO PRO	SP3110
Display of 6 TEMPO PRO	SP3110P

SPENGLER QUALITY

CE MARKING

European directive 93/42/CEE : non-contact thermometer class IIa. Spengler non contact thermometers meet the essential requirements of safety health, hygiene and environmental protection.

After sales service

Spengler provides maintenance and repair of equipment

Legal Guarantee

(art. 1641 and 1648 of the Civil Code): Spengler non-contact thermometers have a lifetime guarantee against any hidden manufacture defect from the purchase date.

Spengler guarantee

Spengler warranty of 2 years (under normal use conditions): free replacement of parts recognised

as defective by standard exchange (function defect from a material, manufacture or design defect) by our customer service department.

Consumable parts (batteries) are not covered by the guarantee. The guarantee does not cover the maintenance defects by the user.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 Aix en Provence | FRANCE

Tel. : **+33 (0)4 42 52 89 80**

Mail : **export@spengler.fr**