

Vaquez-Laubry[®] Classic

The aneroid cuff-mounted shockproof sphygmomanometer



MANOMETER

The tensiometer is made up of an indissociable set of mechanisms. Protected against shocks by a two-storey structure. The pressure gauge is supported by full grain leather.



HIGH PRECISION

Graduations/precision needle: measurements readable from all angles.
Dilated measurement scale = increased accuracy.



PROTECTED FROM SHOCKS

Large dial in non-deforming material, anti-corrosion, impact-resistant (in accordance with EN ISO 81060 chap. 6.4). Optimal protection of the mechanism.
Foam pressure gauge protector.



INFLATION BULB

Spengler's "Classic" model. Maximum volume and ergonomic shape: comfortable handling and rapid inflation.



CUFF

Made of undeformable polycotton washable to 30-40° C. Velcro closure (large adult, adult, child) or triple straps (adult). Calibrated (arm and artery circumference boundary markers).



AIR RELEASE VALVE

Progressive decompression controlled by ultra-precise valve.



INFLATABLE BLADDER

Double tubing. Rubber pocket "high resistance". Adapted to the circumference of the arm to ensure accurate tension measurement.



SPHYGMOMANOMETER

Vaquez-Laubry® Classic

Type	aneroid cuff-mounted	
Manometer	housing	zinc alloy, chrome plated Ø 57 mm
	mechanism	Precision: ± 3 mmHg Copper-beryllium alloy with constant elasticity
	dial	graduations : 0-300 mmHg Ø 50 mm markings: item serial number
	protective bumper	blue or black removable rubber
Manometer holder	fine leather	
Complete bulb	bulb	Spengler "Classic" type black PVC Ø 50 mm
	air release valve	"Classic" type chrome plated linear and accurate deflation
Cuff	fastener	velcro or straps
	fabric	polycotton
	markings	index range (size, artery), washing instructions
Bladder/Tube	double tubing	rubber, highly resistant to ozone and UV
	markings	index range (size, artery), washing instructions
	bladder protection	for straps model
Packaging	carrybag	nylon zipper pouch
Spengler guarantee	manometer : 2 years, sauf consommables (bulbs, bladders, cuffs)	
Legale guarantee	For life against any hidden fault of manufacturing	

References	Model	Bladder size
518002	velcro, large adult, dark blue (L)	16 x 32 cm
518001	velcro, adult, grey (M)	12 x 27 cm
518012	Velcro, adult, dark blue (M)	12 x 27 cm
518013	Velcro, red, adult (M)	12 x 27 cm
518021	velcro, red, adulte (M)	12 x 27 cm
518022	velcro, grey, child/small adult (S)	8 x 23 cm
518112	straps, adult, dark blue (M)	12 x 27 cm
518113	straps, adult, red (M)	12 x 27 cm

Contrôle métrologique every two years (European recommendation)

After-sales service : all spare parts are available from your Spengler dealer.



SPENGLER QUALITY

CE Marking

The Spengler tensiometers, including the Vaquez Laubry® are medical devices of class Im measuring function. They comply with current regulations by perfectly meeting the general safety and performance requirements.

Customer technical support center

SPENGLER MEDICAL provides devices maintenance and repair services.
Contact our customer technical support center on +33 4 42 90 31 31 or click on SAV tab on www.spengler.fr

Spengler warranty

Spengler 2-year warranty (under normal conditions of use): repair or replacement free of charge of defective parts recognized (malfunction resulting from a defect in material, manufacturing or design) by our after-sales service. The armbands are guaranteed for 12 months. The warranty does not cover maintenance defects on the part of the user.



www.spengler.fr

Spengler S.A.S | 30 rue Jean de Guiramand | 13290 AIX-EN-PROVENCE | FRANCE

Tel. : **+33 (0)4 42 90 31 31**
Mail : **contact@spengler.fr**