

## DESJARDIN DUCT DILATOR N°1

Réf. IE40201



## DESJARDIN DUCT DILATOR N°1



Réf. IE40201

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040017648
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,9	Width packaging (cm):	
Length (cm) :	37,4	Length packaging (cm) :	
Height (cm) :	0,4	Height packaging (cm) :	
Weight (kg) :	0,027	Weight packaging (kg) :	0,027
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 21:06 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IE40201

All characteristics are given for information only and may be changed