

## MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 16 CM

Réf. IP49016



# MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 16 CM



Réf. IP49016

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040138305
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	22,5
Height (cm) :		Height packaging (cm) :	1,0
Weight (kg) :	0,070	Weight packaging (kg) :	0,070
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:49 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP49016

All characteristics are given for information only and may be changed