

MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 18 CM

Réf. IP49018



MIXTER-O'SHAUGNESSY ARTERY FORCEPS 90° 18 CM



Réf. IP49018

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040138312
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,5
Length (cm) :		Length packaging (cm) :	19,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :	0,070	Weight packaging (kg) :	0,070
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:49 - This article is available on our website : www.gsh-med.fr with the following reference IP49018

All characteristics are given for information only and may be changed