

SAUERBRÜCK BONE RONGEUR MULTIACTION STRAIGHT CURVED JAW 30 CM

Réf. IP49430



SAUERBRÜCK BONE RONGEUR MULTIACTION STRAIGHT CURVED JAW 30 CM



Réf. IP49430

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040027388
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	17,5	Width packaging (cm):	
Length (cm) :	31,0	Length packaging (cm) :	
Height (cm) :	1,5	Height packaging (cm) :	
Weight (kg) :	0,643	Weight packaging (kg) :	0,643
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:23 - This article is available on our website : www.gsh-med.fr with the following reference IP49430

All characteristics are given for information only and may be changed