

BULLDOG CLAMP STRAIGHT 5 CM

Réf. IP55005



BULLDOG CLAMP STRAIGHT 5 CM



Réf. IP55005

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040028385
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,7	Width packaging (cm):	
Length (cm) :	5,3	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,008	Weight packaging (kg) :	0,008
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:27 - This article is available on our website : www.gsh-med.fr with the following reference IP55005

All characteristics are given for information only and may be changed