

BULLDOG CLAMP DE BAKEY STRAIGHT 5 CM

Réf. IP55105



BULLDOG CLAMP DE BAKEY STRAIGHT 5 CM



Réf. IP55105

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040028422
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,0	Width packaging (cm):	
Length (cm) :	5,0	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,006	Weight packaging (kg) :	0,008
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:27 - This article is available on our website : www.gsh-med.fr with the following reference IP55105

All characteristics are given for information only and may be changed