

## BULLDOG CLAMP DE BAKEY STRAIGHT 6CM

Réf. IP55106



## BULLDOG CLAMP DE BAKEY STRAIGHT 6CM



Réf. IP55106

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040028439
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,4	Width packaging (cm):	
Length (cm) :	16,9	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,011	Weight packaging (kg) :	0,011
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:27 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP55106

All characteristics are given for information only and may be changed