

## BULLDOG CLAMP DE BAKEY STRAIGHT 7.5 CM

Réf. IP55107



# BULLDOG CLAMP DE BAKEY STRAIGHT 7.5 CM



Réf. IP55107

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE	Barcode	3661040028446
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	13,0	Width packaging (cm):	
Length (cm) :	6,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,015	Weight packaging (kg) :	0,015
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:27 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP55107

All characteristics are given for information only and may be changed