

## BULLDOG CLAMP DE BAKEY STRAIGHT 12 CM

Réf. IP55112



# BULLDOG CLAMP DE BAKEY STRAIGHT 12 CM



Réf. IP55112

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE        | Barcode                 | 3661040028477 |
| Guarantee (month):         | 24        | Customs Code            | 9018.9090     |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               | 7,0       | Width packaging (cm):   |               |
| Length (cm) :              | 20,5      | Length packaging (cm) : |               |
| Height (cm) :              | 1,0       | Height packaging (cm) : |               |
| Weight (kg) :              | 0,020     | Weight packaging (kg) : | 0,020         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



19/05/2022 21:16 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP55112

All characteristics are given for information only and may be changed