

## DIALYSYS FORCEPS 16 CM

Réf. IP75416



## DIALYSYS FORCEPS 16 CM



Réf. IP75416

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040028644
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,0	Width packaging (cm):	
Length (cm) :	19,0	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,047	Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:29 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP75416

All characteristics are given for information only and may be changed