

DODEN IMPLANTATION FORCEPS 12.5 CM

Réf. IP52201



DODEN IMPLANTATION FORCEPS 12.5 CM



Réf. IP52201

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040138237
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	6,6	Width packaging (cm):	6,6
Length (cm) :	20,6	Length packaging (cm) :	20,6
Height (cm) :	0,7	Height packaging (cm) :	0,7
Weight (kg) :	0,059	Weight packaging (kg) :	0,059
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:40 - This article is available on our website : www.gsh-med.fr with the following reference IP52201

All characteristics are given for information only and may be changed