

## DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 7 CM STRAIGHT

Réf. IP54007



# DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 7 CM STRAIGHT



Réf. IP54007

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040028200
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :	7,0	Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,017	Weight packaging (kg) :	0,017
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



15/11/2022 12:44 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP54007

All characteristics are given for information only and may be changed