


DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 7 CM CURVED

Réf. IP54017



DE WECKER HEMOSTATIC FORCEPS SERREFINE FOR SUTURE 7 CM CURVED



Réf. IP54017

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040028224
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	0,8	Width packaging (cm):	
Length (cm) :	6,8	Length packaging (cm) :	
Height (cm) :	0,7	Height packaging (cm) :	
Weight (kg) :	0,040	Weight packaging (kg) :	0,040
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 22:04 - This article is available on our website : www.gsh-med.fr with the following reference IP54017

All characteristics are given for information only and may be changed