

## HUDSON BONE HAND DRILL WITH 5 DRILLS AND 1 EXTENSIONS

Réf. IP78100



# HUDSON BONE HAND DRILL WITH 5 DRILLS AND 1 EXTENSIONS



Réf. IP78100

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040028705
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	14,2	Width packaging (cm):	
Length (cm) :	28,3	Length packaging (cm) :	
Height (cm) :	4,4	Height packaging (cm) :	
Weight (kg) :	1,032	Weight packaging (kg) :	1,032
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



11/08/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP78100

All characteristics are given for information only and may be changed