

WILLIAMS FORCEPS 16 CM

Réf. IP46731



WILLIAMS FORCEPS 16 CM



Réf. IP46731

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040086743
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,0
Length (cm) :		Length packaging (cm) :	17,5
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :	0,038	Weight packaging (kg) :	0,038
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:15 - This article is available on our website : www.gsh-med.fr with the following reference IP46731

All characteristics are given for information only and may be changed