

LEXER GOUGE - 15 MM

Réf. IG01130



LEXER GOUGE - 15 MM



Réf. IG01130

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040039251
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	6,0
Length (cm) :		Length packaging (cm) :	27,0
Height (cm) :		Height packaging (cm) :	4,0
Weight (kg) :		Weight packaging (kg) :	0,168
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)

30/03/2022 16:44 - This article is available on our website : www.gsh-med.fr with the following reference IG01130

All characteristics are given for information only and may be changed