

LEXER GOUGE 25CM 6 MM

Réf. IG01132



LEXER GOUGE 25CM 6 MM



Réf. IG01132

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040044019
Guarantee (month):	24	Customs Code	9018908500
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,0	Width packaging (cm):	
Length (cm) :	25,4	Length packaging (cm) :	
Height (cm) :	2,0	Height packaging (cm) :	
Weight (kg) :	0,131	Weight packaging (kg) :	0,131
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 16:44 - This article is available on our website : www.gsh-med.fr with the following reference IG01132

All characteristics are given for information only and may be changed