

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 200 MM - 2 MM

Réf. IP46844



FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 200 MM - 2 MM



Réf. **IP46844**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026442
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	12,3	Width packaging (cm):	12,3
Length (cm) :	29,2	Length packaging (cm) :	29,2
Height (cm) :	0,7	Height packaging (cm) :	0,7
Weight (kg) :	0,207	Weight packaging (kg) :	0,207
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:16 - This article is available on our website : www.gsh-med.fr with the following reference IP46844

All characteristics are given for information only and may be changed