

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 180 MM - 5 MM

Réf. IP46848



FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 130° - 180 MM - 5 MM



Réf. **IP46848**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026480
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	25,0
Height (cm) :		Height packaging (cm) :	1,5
Weight (kg) :		Weight packaging (kg) :	0,239
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	(N/A)
Humidity min		Humidity max	(N/A)



30/03/2022 21:16 - This article is available on our website : www.gsh-med.fr with the following reference IP46848

All characteristics are given for information only and may be changed