

FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 90° - 180 MM - 2 MM

Réf. IP46850



FERRIS-SMITH-KERRISSON FORCEPS UPWARD CUTTING 90° - 180 MM - 2 MM



Réf. **IP46850**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026503
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	18,0	Width packaging (cm):	18,0
Length (cm) :	29,5	Length packaging (cm) :	29,5
Height (cm) :	0,6	Height packaging (cm) :	0,6
Weight (kg) :	0,200	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:16 - This article is available on our website : www.gsh-med.fr with the following reference IP46850

All characteristics are given for information only and may be changed