

SCHAEDEL TOWEL FORCEPS 9 CM

Réf. IP47009



SCHAEDEL TOWEL FORCEPS 9 CM



Réf. **IP47009**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040026749
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	1,2	Width packaging (cm):	
Length (cm) :	9,0	Length packaging (cm) :	
Height (cm) :	0,9	Height packaging (cm) :	
Weight (kg) :	0,016	Weight packaging (kg) :	0,016
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:19 - This article is available on our website : www.gsh-med.fr with the following reference IP47009

All characteristics are given for information only and may be changed