


LUER BONE RONGEUR MULTIACTION STRAIGHT 27 CM

Réf. IP48427



LUER BONE RONGEUR MULTIACTION STRAIGHT 27 CM



Réf. IP48427

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040027111
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	9,6	Width packaging (cm):	
Length (cm) :	17,6	Length packaging (cm) :	
Height (cm) :	0,9	Height packaging (cm) :	
Weight (kg) :	0,413	Weight packaging (kg) :	0,413
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:21 - This article is available on our website : www.gsh-med.fr with the following reference IP48427

All characteristics are given for information only and may be changed