

SPLINT FORCEPS 8 CM WITH MAGNIFIER (DIAM 30 MM)

Réf. IP01400



SPLINT FORCEPS 8 CM WITH MAGNIFIER (DIAM 30 MM)



Réf. IP01400

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040021362
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	3,0	Width packaging (cm):	
Length (cm) :	8,0	Length packaging (cm) :	
Height (cm) :	0,9	Height packaging (cm) :	
Weight (kg) :	0,020	Weight packaging (kg) :	0,020
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:13 - This article is available on our website : www.gsh-med.fr with the following reference IP01400

All characteristics are given for information only and may be changed