

LUC NASAL FORCEPS PERFORATED MEDIUM SIZE 19 CM

Réf. IP19719



LUC NASAL FORCEPS PERFORATED MEDIUM SIZE 19 CM



Réf. IP19719

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040023441
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,5	Width packaging (cm):	
Length (cm) :	23,5	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,060	Weight packaging (kg) :	0,060
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 17:06 - This article is available on our website : www.gsh-med.fr with the following reference IP19719

All characteristics are given for information only and may be changed