

## HUNTER SPLINT FORCEPS STRAIGHT 10.5 CM

Réf. IP45510



# HUNTER SPLINT FORCEPS STRAIGHT 10.5 CM



Réf. IP45510

FEATURES			
Medical device	Yes	(Sans)	
Marquage	AUCUN	Barcode	3661040025803
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	0,012	Weight packaging (kg) :	0,012
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:12 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP45510

All characteristics are given for information only and may be changed