

RAMPLAY SPONGE FORCEPS STRAIGHT 20 CM

Réf. IP45617



RAMPLAY SPONGE FORCEPS STRAIGHT 20 CM



Réf. IP45617

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040025841
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,2	Width packaging (cm):	
Length (cm) :	28,6	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,051	Weight packaging (kg) :	0,051
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



19/05/2022 21:15 - This article is available on our website : www.gsh-med.fr with the following reference IP45617

All characteristics are given for information only and may be changed