

SPURLING LAMINECTOMY RONGEUR FORCEPS STRAIGHT 18 CM 4 X 10 MM

Réf. IP51118



SPURLING LAMINECTOMY RONGEUR FORCEPS STRAIGHT 18 CM 4 X 10 MM



Réf. **IP51118**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040027739
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,5	Width packaging (cm):	
Length (cm) :	25,1	Length packaging (cm) :	
Height (cm) :	0,5	Height packaging (cm) :	
Weight (kg) :	0,080	Weight packaging (kg) :	0,080
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/05/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference IP51118

All characteristics are given for information only and may be changed