

## MASSON NEEDLE HOLDER 27 CM

Réf. IP28427



# MASSON NEEDLE HOLDER 27 CM



Réf. **IP28427**

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe I                |               |
| Marquage                   | CE        | Barcode                 | 3661040024691 |
| Guarantee (month):         | 24        | Customs Code            | 90189090      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              | 0,102     | Weight packaging (kg) : | 0,102         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



24/10/2022 21:08 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference IP28427

All characteristics are given for information only and may be changed