

OLSEN-HÉGAR NEEDLE HOLDER 14 CM

Réf. IP32414



OLSEN-HÉGAR NEEDLE HOLDER 14 CM



Réf. IP32414

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | Yes | Classe Ir | |
| Marquage | CE | Barcode | 3661040024936 |
| Guarantee (month): | 24 | Customs Code | 9018.9090 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | |
| Length (cm) : | 14,0 | Length packaging (cm) : | |
| Height (cm) : | | Height packaging (cm) : | |
| Weight (kg) : | 0,034 | Weight packaging (kg) : | 0,034 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



30/03/2022 17:17 - This article is available on our website : www.gsh-med.fr with the following reference IP32414

All characteristics are given for information only and may be changed