

SCHWARZ BONE RONGEUR SIMPLE ARTICULATION 14.5 CM X 3 MM

Réf. IP48014



SCHWARZ BONE RONGEUR SIMPLE ARTICULATION 14.5 CM X 3 MM



Réf. **IP48014**

FEATURES			
Medical device	Yes	Classe Ir	
Marquage	CE	Barcode	3661040106540
Guarantee (month):	24	Customs Code	9018.9090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	8,0	Width packaging (cm):	
Length (cm) :	18,0	Length packaging (cm) :	
Height (cm) :	1,0	Height packaging (cm) :	
Weight (kg) :	0,133	Weight packaging (kg) :	0,133
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/03/2022 21:21 - This article is available on our website : www.gsh-med.fr with the following reference IP48014

All characteristics are given for information only and may be changed